


Product Catalog

2024

Implant Solutions PREMIUM CAD/CAM PRODUCT



Implant Library System
(3Shape, exocad, Dentalwings)

LL/LS [Multibase]

CMFit [Premill, Blank]

3D-Analog

SC / GMS / GTS [Scanbody]

Sirona Multibase

CCM Casting Abutment

Geo Screw Driver

GeOMEDI

We are exploring excellent dental implant solutions

Since 2009, we have been developing patient-customize abutment, and all of our dentists and dental technicians are also trying to make the fittest implant prosthesis for each individual person.

For Ideal Implant Restorations

We have a commitment to provide the best implant solution at a reasonable price, easy and fast way for everyone.

Co-Cr Abutment Wax up & Casting

CAD/CAM Abutment CAD design & CAM milling

For Digital Implant Solution

It is possible to make implant prosthesis by using a scanner, dental CAD, milling machine and 3D printer as developing our own digital solution for dental clinic and laboratory.

Contents

■ Company History / Certification	04
■ Dental CAD/CAM System Workflow	05
■ Design concept	06
GMS Scanbody /SC Scanbody	06
TiN screw / Geo LS Multibase	08
Geo LL Multibase	09
Geo CMFit / Reverse Jig	10
Geo Multi-Holder for CMFit	11
3D-Analog	12
Geo CAD/CAM prosthetic system	13
with digital library	
■ Compatible Implant Systems	14
Osstem Implant TS	14
Osstem Implant US	15
Osstem Implant SS	16
Dentium Super-Line / Implantium	17
Dentium NR-Line	18
Megagen AnyRidge	19
Megagen AnyOne	20
Megagen EZ-plus	21
Megagen Ex-Feel [Internal]	22
Megagen MiNi	22
DIO SM	23
DIO UF / UFII	23
NeoBiotech IS	24
Dentis S-Clean / OneQ	25
Straumann Bone Level	26
Straumann Tissue Level	27
Straumann BLX	27
Nobel Biocare Active	28
Nobel Biocare Branemark	29
Nobel Biocare Replace	30
Dentsply Xive / Frialit2	31
Dentsply Ankylos	32
Biomet 3i Certain	33
Biomet 3i Osseotite	34
Implant Direct 3.0	35
Zimmer Tapered Screw-vent	36
Zimmer Spline	37
Astra Tech Osseospeed TX	38
Astra Tech Osseospeed EV	39
Camlog	40
Conelog	40
MIS SEVEN	41
MIS C1	42
SIC ace / max	43
Thommen SPI	44
BEGO	45
MYTIS Arrow Implant	46
Kyocera POI	47
Kyocera FINESIA	47
Biohorizons Laser-Lok 3.0	48
Anthogyr AXIOM	49

Warantec Oneplant	50
Warantec IU	50
Tekka In-Kone	51
Keystone PrimaConnex / Max-TL	52
Neodent CM	53
Neodent GM	53
Bicon Implant	54
Medentika IPS	55
GC Aadvia Implant	55
BIOTECH KONTACT	56
Medentis medical ICX	57
B&B Dental	58
Bredent Medical Sky	58
NEOSS ProActive	59
Champions (R)evolution	59

■ MultiUnit(Multi-base)	60
■ GTS Transfer Scanbody	62
■ Sirona Multibase	64
■ CCM Casting Abutment	68
■ Geo Screw Driver	73
■ Instructions of CAD/CAM Abutment	76
■ Instructions of CCM Casting Abutment	77
■ Geo Library	79
■ Selection of Library	80
■ Geo Bridge-Single Library	81
■ Products List	83

Company History & Certification

- 2009 ● 01 Released Geo CCM Abutment
- 2010 ● 03 Established Geodent.
- 2011 ● 12 Released Geo Link, CMFit, Library, Geo Driver
- 2012 ● 10 The 1st Geo Symposium at Seoul National University
- 2013 ● 08 Obtained KGMP
08 Obtained KFDA – Implant, endosseous, laboratory supplies
10 Obtained KFDA – Geo Link abutment, Geo CCM abutment
12 Obtained KFDA – Geo CMFit abutment
12 Established GeoMedi Co.,Ltd.
- 2014 ● 07 Obtained KFDA – TCP
09 The 2nd Geo Symposium at Seoul Trade Exhibition & Convention
- 2015 ● 02 Obtained ISO13485, International Organization for Standardization' Medical devices Certification
02 Obtained CE, the European quality certification
04 Released Geo CAD/CAM system, Geo Mill
06 Certified as a venture business
07 Expansion & Relocation of corporation
09 Established research institute affiliated with a corporation
- 2016 ● 03 Obtained KFDA – Geo Scanbody, Geo 3D Analog, Geo screw
06 Obtained KFDA – Geo Multibase abutment
09 Obtained CFDA (China) – Geo Scanbody, Geo 3D analog
- 2017 ● 07 Obtained CE – Geo CCM abutment, Geo CMFit abutment *(Withdrawal, 2023.05)
07 Obtained EN ISO13485 : 2012
11 Obtained certification of INNOBIZ – Ministry of SMEs and Startups
- 2018 ● 05 Obtained CE Multibase – Geo Multibase abutment
11 Obtained EN ISO13485 : 2016 (31720001)
12 Achieved 1 Million dollors in export
- 2019 ● 01 Obtained FDA – Keystone connection abutment (Link, Premill)
09 Obtained certification of MainBIZ – Ministry of SMEs and Startups
- 2020 ● 02 Obtained KFDA – Geo GTS transfer scanbody, Geo Multiscrew abutment
10 Obtained PMDA (Japan) – Geo Multibase abutment AR(SUM/SURO/ASR/ASW)
- 2021 ● 03 Obtained PMDA (Japan) – Geo Multibase abutment IR(ITR/ITW/ISR/ISW)
04 Being selected as a Small-Giant enterprises by Ministry of Employment and Labor
04 Obtained Kyrgyzstan – Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment, Geo Multiscrew abutment, Geo GTS transfer scanbody, Geo healing cap, Geo Scanbody, Geo Impression coping, Geo 3D analog
05 Obtained KFDA – OMEGA System
05 Obtained CE – Geo GTS transfer scanbody, Geo Multiscrew abutment *(Withdrawal, 2023.05)
12 Certified as a venture business – Innovative Growth
- 2022 ● 04 Obtained RZN (Russia) – Geo GMS scanbody
04 Obtained PMDA (Japan) – Geo Multibase abutment BD(BLS/BLM/BLR)
07 Obtained Kazakhstan – Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment, Geo Multiscrew abutment, Geo healing cap, Geo Scanbody, Geo GMS scanbody, Geo Impression coping, Geo 3D analog, Geo LB Cylinder, Geo Temporary Cylinder
11 Obtained RZN (Russia) – Geo Multiscrew abutment, Geo CCM abutment, Geo CMFit abutment, Geo Multibase abutment
11 Obtained FDA – Biomet 3i OSSEOTITE² Certain², NobelActive, NobelReplace Conical Connection, NobelParallel Conical, Nobel Replace Internal Tri-Channel Connection abutment (Link, Premill)
- 2023 ● 10 Obtained Saudi Arabia – Geo Multibase abutment
- 2024 ● 05 Obtained Thailand – Geo Scanbody, Geo GMS scanbody, Geo Impression coping, Geo 3D analog, Geo screwdriver, Geo Torque wrench, Geo Screw

Dental CAD/CAM System Workflow

1. Scanbody inserting




Model and Intra-oral

2. Scanning



3. Designing - 3shape / exocad



Geo Library System

Geo
Abutment & Libraries
with CAD/CAM

7. Setting



6. Sintering

5. RP and 3D printing



3D Analog

4. Milling



Multi-jig for CMFit



Titanium blank

Design concept

GMS Scanbody (Intra-oral)

- Biocompatible Titanium (Ti-6Al-4V ELI) from CARPENTER Co. guarantees high-durability.
- Autoclavable.
- Special coating allows spray-free scanning.
- Providing two different sizes - 7.5mm(Short) / 11.5mm(Long)
- Consistent performance & Visibility in X-ray.
- Various implant systems are available.
- Providing highly precise Geo Library for various CAD software. (3Shape, Exocad, Dental Wings)
- Customized abutment design for CNC & Pre-milled blank is available.
- Zirconia prosthesis designing with multi-base.



Design concept

GMS Scanbody (Intra-oral)




SC Scanbody(Model)




Design concept

TiN Screw

- Tin-coated (Surface treatment)
- High abrasion resistance
 - Improved fixation of screw



LS Multibase



Design concept

LL Multibase

LL Multibase height is editable on Geo Library.



Tin-coating - High esthetic.
LL Multibase for Posterior.

Grooves - Good adhesive strength.
No need for sandblasting.

Natural curve to make emergence profile



Zirconia - Screw retained restorations
Customized Zirconia abutment
Geo Multibase & Screw made of Titanium Grade ELI
Two Screws - TiN coating screw & Titanium screw



Design concept

Geo CMFit [Premill Blank]



Two types of diameter size are available ($\Phi 10$, $\Phi 14$)

Laser marking - Geo implant system code

Two screws - Tin coating screw & Titanium screw

Qualified material - Ti-6Al-4V ELI made by CARPENTER co.

Reverse Jig



The capability of occlusal milling


No support pin is required.

Original screws are used to fix premill blank.

* GEO milling holder & machine flexibility are required.

Design concept

Geo Multi-Holder for CMFit



● Datron



● imes-icore

Multi-holder and Reverse holder





Geo multi holder shows compatibility with various milling machines.
Available with all geo premill blanks.

● VHF Z4




● VHF N4




Design concept

Geo 3D Analog



Pioneer of providing 3D digital Analog
Can be used both stone model & 3D printed model
Reasonable price with excellent durability
Provided all type of implant system
Stainless Steel



Design concept

Geo CAD/CAM prosthetic system with digital library

Dental CAD s/w

3shape

exocad

dental wings


Dental CAM s/w

FOLLOW-ME!
TECHNOLOGY GROUP
hyperDENT

sum
Geo step beyond

vhf

Geo Scanbody




Scanbody that provides the precise library matching according to the user's environment.

CMFit Abutment

Geo libraries provide stable and accurate sitting of customized abutments



Multibase Abutment




Geo Ti-base abutment is well known for its emphasis on aesthetic and functional design.

*Various gingival heights are available with dedicated library.


3D Analog

Multi-functional analog can be used for both plaster(stone) model and 3D printing model.



Retention groove on 3D analog and extra space to fill up wax provides stable fix

Bridge & Single Library




By cooperation with dental CAM s/w companies, precise direct interface milling from Geo-geometry is available.

Compatible Implant Systems

Osstem Implant TS

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase				
	 (11.5mm)  (7.5mm)							
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)				
				 1.2Hex <table><tr><td>20Ncm</td><td>Ø3.0 / 3.5</td></tr><tr><td>30Ncm</td><td>Ø4.0/4.5/5.0</td></tr></table>	20Ncm	Ø3.0 / 3.5	30Ncm	Ø4.0/4.5/5.0
20Ncm	Ø3.0 / 3.5							
30Ncm	Ø4.0/4.5/5.0							

Description-Library code		0101 - SUM	0102 - SURO
Product	Connection Platform Ø	TS Mini  Ø3.0 / 3.5	TS Standard  Ø4.0/4.5/5.0
SC-Scanbody		SC-SUM	SC-SURO
GMS-Scanbody(Short)		GMS-SUMS	GMS-SUROS
GMS-Scanbody(Long)		GMS-SUML	GMS-SUROL
LL-Multibase (Single)	GH : 1.3	LL2-SUM13-H2	LL2-SURO13-H2
	GH : 2.0	LL2-SUM20-H2	LL2-SURO20-H2
	GH : 3.0	-	LL2-SURO30-H2
	GH : 4.0	-	LL2-SURO40-H2
LL-Multibase (Bridge)	GH : 1.3	LL2-SUM13-N2	LL2-SURO13-N2
	GH : 2.0	LL2-SUM20-N2	LL2-SURO20-N2
	GH : 3.0	-	LL2-SURO30-N2
	GH : 4.0	-	LL2-SURO40-N2
LS-Multibase / Single		LS2-SUM-H2	LS2-SURO-H2
LS-Multibase / Bridge		LS2-SUM-N2	LS2-SURO-N2
CMFit Abutment(S) / Ø10		CMFit10-SUM-H2	CMFit10-SURO-H2
CMFit Abutment(B) / Ø10		CMFit10-SUM-N2	CMFit10-SURO-N2
CMFit Abutment(S) / Ø14		CMFit14-SUM-H2	CMFit14-SURO-H2
CMFit Abutment(B) / Ø14		-	CMFit14-SURO-N2
3D-Analog		3D-SUM	3D-SURO
Screw(5ea)		AC-SUM5	AC-SURO5

Compatible Implant Systems

Osstem Implant US

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				
				1.2Hex
				30Ncm


Osstem Implant US

Description-Library code	0103 - EXM	0104 - EXR	0105 - EXW
Connection Platform Ø	US Mini  Ø3.5	US Regular  Ø4.1	US Wide  Ø5.1
Product			
SC-Scanbody	SC-EXM	SC-EXR	SC-EXW
GMS-Scanbody(Short)	GMS-EXMS	GMS-EXRS	GMS-EXWS
GMS-Scanbody(Long)	GMS-EXML	GMS-EXRL	GMS-EXWL
LL-Multibase / Single	LL2-EXM-H2	LL2-EXR-H2	LL2-EXW-H2
LL-Multibase / Bridge	LL2-EXM-N2	LL2-EXR-N2	LL2-EXW-N2
LS-Multibase / Single	LS2-EXM-H2	LS2-EXR-H2	LS2-EXW-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXM-H2	CMFit10-EXR-H2	CMFit10-EXW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-
CMFit Abutment(S) / Ø14	CMFit14-EXM-H2	CMFit14-EXR-H2	CMFit14-EXW-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(Sea)	AC-EXM5	AC-EXR5	AC-EXW5

Compatible Implant Systems

Osstem Implant SS

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (10mm) (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex 30Ncm

Description-Library code	0106 - ISR	0107 - ISW
Connection Platform Ø	SS Regular  Ø4.8	SS Wide  Ø6.0
Product		
SC-Scanbody	SC-ISRW	SC-ISRW
GMS-Scanbody(Short)	GMS-ISRS	GMS-ISWS
GMS-Scanbody(Long)	GMS-ISRL	GMS-ISWL
LL-Multibase / Single	LL2-ISR-02	LL2-ISW-02
LL-Multibase / Bridge	LL2-ISR-N2	LL2-ISW-N2
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-ISR-02	CMFit10-ISW-02
CMFit Abutment / Ø14	CMFit14-ISR-02	CMFit14-ISW-02
3D-Analog	3D-ISR	3D-ISW
Screw(5ea)	AC-EXR5	AC-EXR5

Compatible Implant Systems

Dentium Super-Line / Implantium

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex 25-30Ncm


Dentium Super-Line / Implantium

Description-Library code		0201 - DER	
Product	Connection Platform Ø	Universal	 Ø3.7/4.0/4.5
SC-Scanbody		SC-DER	
GMS-Scanbody(Short)		GMS-DERS	
GMS-Scanbody(Long)		GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DER13-H2	
	GH : 2.0	LL2-DER20-H2	
	GH : 3.0	LL2-DER30-H2	
LL-Multibase (Bridge)	GH : 1.3	LL2-DER13-N2	
	GH : 2.0	LL2-DER20-N2	
	GH : 3.0	LL2-DER30-N2	
LS-Multibase / Single		LS2-DER-H2	
LS-Multibase / Bridge		LS2-DER-N2	
CMFit Abutment(S) / Ø10		CMFit10-DER-H2	
CMFit Abutment(B) / Ø10		CMFit10-DER-N2	
CMFit Abutment(S) / Ø14		CMFit14-DER-H2	
CMFit Abutment(B) / Ø14		CMFit14-DER-N2	
3D-Analog		3D-DER	
Screw(5ea)		AC-DER5	

Compatible Implant Systems

Dentium NR-Line

*Set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
		-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-	-		 1.0square 20Ncm

Description—Library code	0202 - DNR		
Product	Connection Platform Ø	Standard	 Ø3.1/3.6/4.3
OS-Scanbody	OS-DNR		
LL-Multibase / Single	-		
LL-Multibase / Bridge	-		
LS-Multibase / Single	-		
CMFit Abutment / Ø10	CMFit10-DNR-R		
CMFit Abutment / Ø14	-		
3D-Analog	-		
Screw(5ea)	AC-DNR5		

Compatible Implant Systems

Megagen AnyRidge

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex 35Ncm

Megagen AnyRidge

Description-Library code		0301 - MARR	
Product	Connection Platform Ø	Universal	 Ø3.5/4.0/4.5
SC-Scanbody		SC-MARR	
GMS-Scanbody(Short)		GMS-MARRS	
GMS-Scanbody(Long)		GMS-MARRL	
LL-Multibase (Single)	GH : 1.3	LL2-MARR-H2	
	GH : 2.0	LL2-MARR20-H2	
	GH : 3.0	LL2-MARR30-H2	
LL-Multibase (Bridge)	GH : 1.3	LL2-MARR-N2	
	GH : 2.0	LL2-MARR20-N2	
	GH : 3.0	-	
LS-Multibase / Single		LS2-MARR-H2	
CMFit Abutment / Ø10		CMFit10-MARR-H2	
CMFit Abutment / Ø14		CMFit14-MARR-H2	
3D-Analog		3D-MARR	
Screw(5ea)		AC-MARR5	

Compatible Implant Systems

Megagen AnyOne

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral) *Using 1.27Hex Driver*	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex 25Ncm

Description—Library code		0302 - DERA	
Product	Connection Platform Ø	Universal	 Ø 3.5/4.0/4.5/5.0
SC-Scanbody		SC-DER	
GMS-Scanbody(Short)		GMS-DERS	
GMS-Scanbody(Long)		GMS-DERL	
LL-Multibase (Single)	GH : 1.3	LL2-DERA13-H2	
	GH : 2.0	LL2-DERA20-H2	
	GH : 3.0	LL2-DERA30-H2	
LL-Multibase (Bridge)	GH : 1.3	LL2-DERA13-N2	
	GH : 2.0	LL2-DERA20-N2	
	GH : 3.0	LL2-DERA30-N2	
LS-Multibase / Single		LS2-DERA-H2	
LS-Multibase / Bridge		-	
CMFit Abutment(S) / Ø10		CMFit10-DERA-H2	
CMFit Abutment(B) / Ø10		CMFit10-DERA-N2	
CMFit Abutment(S) / Ø14		CMFit14-DERA-H2	
CMFit Abutment(B) / Ø14		CMFit14-DERA-N2	
3D-Analog		3D-DER	
Screw(5ea)		AC-MA05	

Compatible Implant Systems

Megagen EZ-plus

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-	-		 1.2Hex
				35Ncm


Megagen EZ-plus

Description-Library code	0303 - MEPM	0304 - MEPR
Connection Platform Ø	 Ø 3.3	 Ø 4.0/4.5/5.0
Product		
SC-Scanbody	SC-MEPM	SC-MEPR
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	-	-
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-MEPM-E	CMFit10-MEPR-E
CMFit Abutment / Ø14	-	-
3D-Analog	-	-
Screw(5ea)	AC-MEPM5	AC-MEPR5

Compatible Implant Systems

Megagen ExFeel Internal


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				30Ncm

Description-Library code	0305 - ITRK
Connection Platform Ø	Regular  Ø4.8
Product	
SC-Scanbody	SC-ITRK
GMS-Scanbody(Short)	-
GMS-Scanbody(Long)	-
LL-Multibase / Single	LL2-ITRK-02
LL-Multibase / Bridge	LL2-ITRK-N2
LS-Multibase / Single	-
CMFit Abutment(S) / Ø10	CMFit10-ITRK-02
CMFit Abutment(B) / Ø10	CMFit10-ITRK-N2
CMFit Abutment(S) / Ø14	CMFit14-ITRK-02
3D-Analog	3D-ITR
Screw(Sea)	AC-EXR5

Megagen MiNi

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Multibase	Ø10 CMFit Abutment	Screw	Driver (Ncm)
	-			 1.2Hex
				15Ncm Mini 3.0/3.25

Description-Library code	0306 - MMN
Connection Platform Ø	MINI  Ø3.0/3.25
Product	
SC-Scanbody	SC-MMN
LL-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-MMN-H
Screw(Sea)	AC-MMN5

Compatible Implant Systems

DIO SM

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	Ø10 CMFit Abutment	Single LS-Multibase	Screw	Driver (Ncm)
		-		 1.7 Torx 30Ncm

Description—Library code	0401 - DSMM	0402 - DSMR
Connection Platform Ø	Mini  Ø 3.8 / 4.1	Regular  Ø 4.5 / 5.0
Product		
Scanbody	SC-DSMM	SC-DSMR
LS-Multibase / Single	-	-
CMFit Abutment Ø10	CMFit10-DSMM-T	CMFit10-DSMR-T
Screw(5ea)	AC-DSMM5	AC-DSMR5

DIO UF / UFII

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm)  (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
			 	 1.2Hex 20Ncm Ø 3.0/3.3 30Ncm Ø 3.8/4.0/4.5/5.0

Description—Library code	0403 - DN30	0404 - SURD
Connection Platform Ø	UFII Narrow  Ø 3.0/3.3	Standard  Ø 3.8/4.0/4.5/5.0
Product		
SC-Scanbody	SC-MMN	SC-SURO
GMS-Scanbody(Short)	-	GMS-SURO5
GMS-Scanbody(Long)	-	GMS-SURO1
LL-Multibase (Single)	GH : 1.3	-
	GH : 2.0	-
	GH : 3.0	-
LL-Multibase (Bridge)	GH : 1.3	-
	GH : 2.0	-
	GH : 3.0	-
CMFit Abutment(S) / Ø10	CMFit10-DN30-H	CMFit10-SURD-H2
CMFit Abutment(B) / Ø10	-	CMFit10-SURD-N2
CMFit Abutment(S) / Ø14	-	CMFit14-SURD-H2
CMFit Abutment(B) / Ø14	-	CMFit14-SURD-N2
3D-Analog	-	3D-SURO
Screw(5ea)	AC-DN305	AC-SURD5

DIO SM

DIO UF / UFII

Compatible Implant Systems

NeoBiotech IS

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				30Ncm Ø 3.5/4.0/4.5/5.0

Description-Library code		0501 - DERN
Product	Connection Platform Ø	Standard  Ø 3.5/4.0/4.5/5.0
SC-Scanbody		SC-SURO
GMS-Scanbody(Short)		GMS-DERNS
GMS-Scanbody(Long)		GMS-DERNL
LL-Multibase (Single)	GH : 1.3	LL2-DERN-H2
	GH : 2.0	LL2-DERN20-H2
LL-Multibase (Single)	GH : 1.3	LL2-DERN-N2
	GH : 2.0	LL2-DERN20-N2
LS-Multibase / Single		-
CMFit Abutment(S) / Ø10		CMFit10-DERN-H2
CMFit Abutment(B) / Ø10		CMFit10-DERN-N2
CMFit Abutment(S) / Ø14		CMFit14-DERN-H2
CMFit Abutment(B) / Ø14		CMFit14-DERN-N2
3D-Analog		3D-DER
Screw(Sea)		AC-DERN5

Compatible Implant Systems

Dentis S-Clean / One Q

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex 30Ncm

Dentis S-Clean / One Q

Description-Library code		0601- DERT
Product	Connection Platform Ø	Regular  Ø 4.2/4.7/5.2/6.0
SC-Scanbody		SC-DER
GMS-Scanbody(Short)		GMS-DERS
GMS-Scanbody(Long)		GMS-DERL
LL-Multibase (Single)	GH : 1.3	LL2-DERT13-H2
	GH : 2.0	-
LL-Multibase (Bridge)	GH : 1.3	LL2-DERT13-N2
	GH : 2.0	-
LS-Multibase / Single		-
CMFit Abutment(S) / Ø10		CMFit10-DERT-H2
CMFit Abutment(S) / Ø14		CMFit14-DERT-H2
3D-Analog		3D-DER
Screw(5ea)		AC-DERT5

Compatible Implant Systems

Straumann Bone Level

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx 35Ncm

Description—Library code		0700 - BLS	0701 - BLM	0702 - BLR
Product	Connection Platform Ø	SC  Ø2.9	NC  Ø3.3	RC  Ø4.1/4.8
SC-Scanbody		-	SC-BLM	SC-BLR
GMS-Scanbody(Short)		-	GMS-BLMS	GMS-BLRS
GMS-Scanbody(Long)		GMS-BLSL	GMS-BLML	GMS-BLRL
LL-Multibase (Single)	GH : 1.3	LL2-BLS-R2	LL2-BLM-R2	LL2-BLR-R2
	GH : 2.0	-	LL2-BLM20-R2	LL2-BLR20-R2
	GH : 3.0	-	-	LL2-BLR30-R2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-BLM-N2	LL2-BLR-N2
	GH : 2.0	-	LL2-BLM20-N2	LL2-BLR20-N2
	GH : 3.0	-	-	-
LS-Multibase / Single		-	LS2-BLM-R2	LS2-BLR-R2
CMFit Abutment / Ø10		CMFit10-BLS-R2	CMFit10-BLM-R2	CMFit10-BLR-R2
CMFit Abutment / Ø14		-	CMFit14-BLM-R2	CMFit14-BLR-R2
3D-Analog		3D-BLS	3D-BLM	3D-BLR
Screw(5ea)		AC-BLS5	AC-BLM5	AC-BLR5

Compatible Implant Systems

Straumann Tissue Level

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				
				Torx
				35Ncm

Description-Library code	0703 - IT35	0704 - ITR	0705 - ITW
Connection Platform Ø	NN  Ø3.5	RN  Ø4.8	WN  Ø6.5
Product			
SC-Scanbody	SC-IT35	SC-ITRW	SC-ITRW
GMS-Scanbody(Short)	-	GMS-ITRS	GMS-ITWS
GMS-Scanbody(Long)	-	GMS-ITRL	GMS-ITWL
LL-Multibase (Single)	GH : 1.3	LL2-ITR-02	LL2-ITW-02
	GH : 2.0	LL2-ITR20-02	LL2-ITW20-02
LL-Multibase (Bridge)	GH : 1.3	LL2-ITR-N2	LL2-ITW-N2
	GH : 2.0	LL2-ITR20-N2	LL2-ITW20-N2
LS-Multibase / Single	LS2-IT35-0	LS2-ITR-02	LS2-ITW-02
CMFit Abutment / Ø10	CMFit10-IT35-0	CMFit10-ITR-02	CMFit10-ITW-02
CMFit Abutment / Ø14	-	CMFit14-ITW-02	CMFit14-ITW-02
3D-Analog	-	3D-ITR	3D-ITW
Screw(5ea)	AC-IT355	AC-IT5	AC-IT5

Straumann BLX

*Set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Single LL-Multibase	∅10 CMFit Abutment	∅14 CMFit Abutment	3D-Analog	Screw																														
 (10mm)																																			
<table><tr><th>Description-Library code</th><td>0710 - BLX</td><td>0711 - BLXW</td></tr><tr><th>Connection Platform ∅</th><td> ∅ 3.75/4.5</td><td> ∅ 5.5/6.5</td></tr><tr><th>Product</th><td></td><td></td></tr><tr><th>GMS-Scanbody(Long)</th><td>GMS-BLXL</td><td>GMS-BLXL</td></tr><tr><th>LL-Multibase / Single</th><td>LL2-BLX-T2</td><td>-</td></tr><tr><th>CMFit Abutment / ∅10</th><td>CMFit10-BLX-T2</td><td>CMFit10-BLXW-T2</td></tr><tr><th>CMFit Abutment / ∅14</th><td>CMFit14-BLX-T2</td><td>CMFit14-BLXW-T2</td></tr><tr><th>3D-Analog</th><td>3D-BLX</td><td>-</td></tr><tr><th>Screw(5ea)</th><td>AC-BLX5</td><td>AC-BLX5</td></tr></table>					Description-Library code	0710 - BLX	0711 - BLXW	Connection Platform ∅	∅ 3.75/4.5	∅ 5.5/6.5	Product			GMS-Scanbody(Long)	GMS-BLXL	GMS-BLXL	LL-Multibase / Single	LL2-BLX-T2	-	CMFit Abutment / ∅10	CMFit10-BLX-T2	CMFit10-BLXW-T2	CMFit Abutment / ∅14	CMFit14-BLX-T2	CMFit14-BLXW-T2	3D-Analog	3D-BLX	-	Screw(5ea)	AC-BLX5	AC-BLX5	<table><tr><th>Driver (Ncm)</th></tr><tr><td> Torx</td></tr><tr><td>35Ncm</td></tr></table>	Driver (Ncm)	 Torx	35Ncm
Description-Library code	0710 - BLX	0711 - BLXW																																	
Connection Platform ∅	∅ 3.75/4.5	∅ 5.5/6.5																																	
Product																																			
GMS-Scanbody(Long)	GMS-BLXL	GMS-BLXL																																	
LL-Multibase / Single	LL2-BLX-T2	-																																	
CMFit Abutment / ∅10	CMFit10-BLX-T2	CMFit10-BLXW-T2																																	
CMFit Abutment / ∅14	CMFit14-BLX-T2	CMFit14-BLXW-T2																																	
3D-Analog	3D-BLX	-																																	
Screw(5ea)	AC-BLX5	AC-BLX5																																	
Driver (Ncm)																																			
 Torx																																			
35Ncm																																			

Compatible Implant Systems

Nobel Biocare Active

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 15Ncm Ø3.0 35Ncm Ø3.5/4.3/5.0/6.0

Description—Library code		0801- AC30	0802- ACN	0803- ACR	0804- ACW
Product	Connection Platform Ø	 Ø3.0	NP  Ø3.5	RP  Ø4.3/5.0	WP  Ø5.5
Scanbody		OS-AC30	SC-ACN	SC-ACR	SC-ACW
GMS-Scanbody(Short)		-	GMS-ACNS	GMS-ACRS	GMS-ACWS
GMS-Scanbody(Long)		-	GMS-ACNL	GMS-ACRL	GMS-ACWL
LL-Multibase (Single)	GH : 1.3	LL2-AC30-T2	LL2-ACN-H2	LL2-ACR-H2	LL2-ACW-H2
	GH : 2.0	-	LL2-ACN20-H2	LL2-ACR20-H2	-
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ACN-N2	LL2-ACR-N2	LL2-ACW-N2
	GH : 2.0	-	LL2-ACN20-N2	LL2-ACR20-N2	-
LS-Multibase / Single		LS2-AC30-T2	LS2-ACN-H2	LS2-ACR-H2	-
CMFit Abutment / Ø10		CMFit10-AC30-T2	CMFit10-ACN-H2	CMFit10-ACR-H2	CMFit10-ACW-H2
CMFit Abutment / Ø14		-	CMFit14-ACN-H2	CMFit14-ACR-H2	CMFit14-ACW-H2
3D-Analog		3D-AC30	3D-ACN	3D-ACR	3D-ACW
Screw(Sea)		AC-AC305	AC-ACN5	AC-ACR5	AC-ACR5

Compatible Implant Systems

Nobel Biocare Branemark

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				35Ncm

Nobel Biocare Branemark

Description-Library code	0805 - EXMB	0806 - EXRB	0807 - EXWB
Connection Platform Ø	NP  Ø3.5	RP  Ø4.1	WP  Ø5.1
Product			
SC-Scanbody	SC-EXMB	SC-EXRB	SC-EXWB
GMS-Scanbody(Short)	-	-	-
GMS-Scanbody(Long)	-	-	-
LL-Multibase / Single	LL2-EXMB-H2	LL2-EXRB-H2	LL2-EXWB-H2
LL-Multibase / Bridge	LL2-EXMB-N2	LL2-EXRB-N2	LL2-EXWB-N2
LS-Multibase / Single	LS2-EXMB-H2	LS2-EXRB-H2	LS2-EXWB-H2
LS-Multibase / Bridge	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EXMB-H2	CMFit10-EXRB-H2	CMFit10-EXWB-H2
CMFit Abutment(S) / Ø14	CMFit14-EXMB-H2	CMFit14-EXRB-H2	CMFit14-EXWB-H2
3D-Analog	3D-EXM	3D-EXR	3D-EXW
Screw(5ea)	AC-EXMB5	AC-EXRB5	AC-EXWB5

Compatible Implant Systems

Nobel Biocare Replace

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip 35Ncm

Description-Library code	0808- RPM	0809- RPR	0810- RPW	0811- RP60
Connection Platform Ø	NP  Ø3.5	RP  Ø4.3	WP  Ø5.0	 Ø6.0
Product				
Scanbody	SC-RPM	SC-RPR	SC-RPW	SC-RP60
GMS-Scanbody(Short)	GMS-RPMS	GMS-RPRS	GMS-RPWS	GMS-RP60S
GMS-Scanbody(Long)	GMS-RPML	GMS-RPRL	GMS-RPWL	GMS-RP60L
LL-Multibase / Single	LL2-RPM-E2	LL2-RPR-E2	LL2-RPW-E2	LL2-RP60-E2
LL-Multibase / Bridge	LL2-RPM-N2	LL2-RPR-N2	LL2-RPW-N2	LL2-RP60-N2
LS-Multibase / Single	LS2-RPM-E2	LS2-RPR-E2	LS2-RPW-E2	LS2-RP60-E2
CMFit Abutment / Ø10	CMFit10-RPM-E2	CMFit10-RPR-E2	CMFit10-RPW-E2	CMFit10-RP60-E2
CMFit Abutment / Ø14	CMFit14-RPM-E2	CMFit14-RPR-E2	CMFit14-RPW-E2	CMFit14-RP60-E2
3D-Analog	3D-RPM	3D-RPR	3D-RPW	3D-RP60
Screw(Sea)	AC-RPM5	AC-RPR5	AC-RPR5	AC-RPR5

Compatible Implant Systems

Dentsply Xive/Frialit2

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.22Hex 24Ncm

Description-Library code	0901 - XV3	0902 - XVM	0903 - XVR	0904 - XVW	0905 - XV55
Connection Platform Ø	 Ø3.0	 Ø3.4	 Ø3.8	 Ø4.5	 Ø5.5
Product					
SC-Scanbody	SC-XV3	SC-XVM	SC-XVR	SC-XVW	SC-XVW
GMS-Scanbody(Short)	-	GMS-XVMS	GMS-XVRS	GMS-XVWS	GMS-XVWS
GMS-Scanbody(Long)	GMS-XV3L	GMS-XVML	GMS-XVRL	GMS-XVWL	GMS-XVWL
LL-Multibase / Single	LL2-XV3-H2	LL2-XVM-H2	LL2-XVR-H2	LL2-XVW-H2	LL2-XV55-H2
LL-Multibase / Bridge	LL2-XV3-N2	LL2-XVM-N2	LL2-XVR-N2	LL2-XVW-N2	LL2-XV55-N2
LS-Multibase / Single	LS2-XV3-H2	LS2-XVM-H2	LS2-XVR-H2	LS2-XVW-H2	LS2-XV55-H2
CMFit Abutment / Ø10	CMFit10-XV3-H2	CMFit10-XVM-H2	CMFit10-XVR-H2	CMFit10-XVW-H2	CMFit10-XV55-H2
CMFit Abutment / Ø14	-	CMFit14-XVM-H2	CMFit14-XVR-H2	CMFit14-XVW-H2	CMFit14-XV55-H2
3D-Analog	3D-XV3	3D-XVM	3D-XVR	3D-XVW	3D-XV55
Screw(Sea)	AC-XV35	AC-XV5	AC-XV5	AC-XV5	AC-XV5

Compatible Implant Systems

Dentsply Ankylos

*Set (Scanbody/Multibase/CMFit Abutment + screw)


OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				1.0Hex 15Ncm

Description—Library code		0906 - AKLL
Product	Connection Platform Ø	universal  Ø 3.4/4.5/5.5
OS-Scanbody		OS-AKLL
LL-Multibase (Single)	GH : 1.3	LL2-AKLL-T
	GH : 2.0	LL2-AKLL20-T
LL-Multibase (Bridge)	GH : 1.3	LL2-AKLL-N
	GH : 2.0	LL2-AKLL20-N
LS-Multibase / Single		-
CMFit Abutment(S) / Ø10		CMFit10-AKLL-T
CMFit Abutment(B) / Ø10		CMFit10-AKLL-N
CMFit Abutment(S) / Ø14		CMFit14-AKLL-T
3D-Analog		3D-AKLL
Screw		-

Compatible Implant Systems

Biomet 3i Certain

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex 20Ncm

Description-Library code	1001 - 3ICM	1002 - 3ICR	1003 - 3ICW	1004 - 3IC60
Connection Platform Ø	 Ø3.4	 Ø4.1	 Ø5.0	 Ø6.0
Product				
SC-Scanbody	SC-3ICM	SC-3ICR	SC-3ICR	SC-3ICR
GMS-Scanbody(Short)	GMS-3ICMS	GMS-3ICRS	GMS-3ICRS	GMS-3ICRS
GMS-Scanbody(Long)	GMS-3ICML	GMS-3ICRL	GMS-3ICRL	GMS-3ICRL
LL-Multibase / Single	LL2-3ICM-H2	LL2-3ICR-H2	LL2-3ICW-H2	LL2-3IC60-H2
LL-Multibase / Bridge	LL2-3ICM-N2	LL2-3ICR-N2	LL2-3ICW-N2	LL2-3IC60-N2
LS-Multibase / Single	LS2-3ICM-H2	LS2-3ICR-H2	LS2-3ICW-H2	LS2-3IC60-H2
CMFit Abutment / Ø10	CMFit10-3ICM-H2	CMFit10-3ICR-H2	CMFit10-3ICW-H2	CMFit10-3IC60-H2
CMFit Abutment / Ø14	-	CMFit14-3ICR-H2	CMFit14-3ICW-H2	CMFit14-3IC60-H2
3D-Analog	3D-3ICM	3D-3ICR	3D-3ICW	3D-3IC60
Screw(Sea)	AC-3IC5	AC-3IC5	AC-3IC5	AC-3IC5

Compatible Implant Systems

Biomet 3i Osseotite

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				
				20Ncm

Description—Library code	1005- EX3M	1006- EXR	1007- EXT	1008- EXTW
Connection Platform Ø	 Ø3.4	 Ø4.1	 Ø5.0	 Ø6.0
Product				
Scanbody	SC-EX3M	SC-EXR	SC-EXR	SC-EXR
GMS-Scanbody(Short)	-	GMS-EXRS	GMS-EXRS	GMS-EXRS
GMS-Scanbody(Long)	-	GMS-EXRL	GMS-EXRL	GMS-EXRL
LL-Multibase / Single	LL2-EX3M-H2	LL2-EXR-H2	LL2-EXT-H2	LL2-EXTW-H2
LL-Multibase / Bridge	LL2-EX3M-N2	LL2-EXR-N2	LL2-EXT-N2	LL2-EXTW-N2
LS-Multibase / Single	-	LS2-EXR-H2	LS2-EXT-H2	-
LS-Multibase / Bridge	-	-	-	-
CMFit Abutment(S) / Ø10	CMFit10-EX3M-H2	CMFit10-EXR-H2	CMFit10-EXT-H2	CMFit10-EXTW-H2
CMFit Abutment(B) / Ø10	-	CMFit10-EXR-N2	-	-
CMFit Abutment(S) / Ø14	-	CMFit14-EXR-H2	CMFit14-EXT-H2	CMFit14-EXTW-H2
3D-Analog	-	3D-EXR	-	-
Screw(Sea)	AC-EXR5	AC-EXR5	AC-EXR5	AC-EXR5

Compatible Implant Systems

Implant Direct 3.0

*Set (Scanbody/Multibase/CMFit Abutment + screw)

OS-Scanbody(Intra-Oral)		Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
			-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.25Hex
				30Ncm

Description—Library code	1101 - TSVM
Connection Platform Ø	Implant Direct  Ø3.0
Product	
OS-Scanbody	OS-TSVM
LL-Multibase / Single	LL2-TSVM-H2
LL-Multibase / Bridge	-
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-TSVM-H2
CMFit Abutment / Ø14	-
3D-Analog	3D-TSVM
Screw(5ea)	AC-TSVM5

Compatible Implant Systems

Zimmer Tapered Screw-vent

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex 30Ncm

Description-Library code		1201- TSV3	1202- TSV4	1203- TSV5
Product	Connection Platform Ø	 Ø3.5	 Ø4.5	 Ø5.7
Scanbody		SC-TSV3	SC-TSV4	SC-TSV5
GMS-Scanbody(Short)		GMS-TSV3S	GMS-TSV4S	GMS-TSV5S
GMS-Scanbody(Long)		GMS-TSV3L	GMS-TSV4L	GMS-TSV5L
LL-Multibase (Single)	GH : 1.0	LL2-TSV3-H2	LL2-TSV4-H2	LL2-TSV5-H2
	GH : 2.0	LL2-TSV320-H2	LL2-TSV420-H2	-
LL-Multibase (Bridge)	GH : 1.0	LL2-TSV3-N2	LL2-TSV4-N2	LL2-TSV5-N2
	GH : 2.0	LL2-TSV320-N2	LL2-TSV420-N2	-
LS-Multibase / Single		LS2-TSV3-H2	LS2-TSV4-H2	LS2-TSV5-H2
CMFit Abutment / Ø10		CMFit10-TSV3-H2	CMFit10-TSV4-H2	CMFit10-TSV5-H2
CMFit Abutment / Ø14		CMFit14-TSV3-H2	CMFit14-TSV4-H2	CMFit14-TSV5-H2
3D-Analog		3D-TSV3	3D-TSV4	3D-TSV5
Screw(5ea)		AC-TSV5	AC-TSV5	AC-TSV5

Compatible Implant Systems

Zimmer Spline

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-		-		 1.25Hex
				30Ncm

Description-Library code	1204 - SP40	1205- SP50
Connection Platform Ø	 Ø 3.7 / 4.0	 Ø 5.0
Product		
SC-Scanbody	SC-SP40	SC-SP40
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SP40-H	LL2-SP50-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	-	-
CMFit Abutment / Ø14	CMFit14-SP40-H	CMFit14-SP50-H
3D-Analog	-	-
Screw(Sea)	AC-TSV5	AC-TSV5

Compatible Implant Systems

Astra Tech Osseospeed TX

*Set (Scanbody/Multibase/CMFit Abutment + screw)



SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase						
	 (11.5mm) (7.5mm)									
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)						
				 1.27Hex <table><tr><td>15Ncm</td><td>Ø3.0</td></tr><tr><td>20Ncm</td><td>Ø3.5/4.0</td></tr><tr><td>25Ncm</td><td>Ø4.5/5.0</td></tr></table>	15Ncm	Ø3.0	20Ncm	Ø3.5/4.0	25Ncm	Ø4.5/5.0
15Ncm	Ø3.0									
20Ncm	Ø3.5/4.0									
25Ncm	Ø4.5/5.0									

Description—Library code		1301- AS3	1302- ASR	1303- ASW
Product	Connection Platform Ø	Yellow  Ø3.0	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0
Scanbody		SC-AS3	SC-ASR	SC-ASW
GMS-Scanbody(Short)		-	GMS-ASRS	GMS-ASWS
GMS-Scanbody(Long)		GMS-AS3L	GMS-ASRL	GMS-ASWL
LL-Multibase (Single)	GH : 1.3	LL2-AS3-T2	LL2-ASR-H2	LL2-ASW-H2
	GH : 2.0	-	LL2-ASR20-H2	LL2-ASW20-H2
LL-Multibase (Bridge)	GH : 1.3	-	LL2-ASR-N2	LL2-ASW-N2
	GH : 2.0	-	LL2-ASR20-N2	LL2-ASW20-N2
LS-Multibase / Single		-	LS2-ASR-H2	LS2-ASW-H2
CMFit Abutment / Ø10		CMFit10-AS3-T2	CMFit10-ASR-H2	CMFit10-ASW-H2
CMFit Abutment / Ø14		-	CMFit14-ASR-H2	CMFit14-ASW-H2
3D-Analog		3D-AS3	3D-ASR	3D-ASW
Screw(5ea)		AC-AS35	AC-ASR5	AC-ASW5

Compatible Implant Systems

Astra Tech Osseospeed EV

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex 25Ncm

Description—Library code		1304 - AEV30	1305 - AEV36	1306 - AEV42	1307 - AEV48	1308 - AEV54
Product	Connection Platform Ø	Green  Ø 3.0	Purple  Ø 3.6	Yellow  Ø 4.2	Blue  Ø 4.8	Brown  Ø 5.4
SC-Scanbody		SC-AEV30	SC-AEV36	SC-AEV42	SC-AEV48	SC-AEV54
GMS-Scanbody(Short)		GMS-AEV30S	GMS-AEV36S	GMS-AEV42S	GMS-AEV48S	GMS-AEV54S
GMS-Scanbody(Long)		GMS-AEV30L	GMS-AEV36L	GMS-AEV42L	GMS-AEV48L	GMS-AEV54L
LL-Multibase (Single)	GH : 1.3	LL2-AEV30-T2	LL2-AEV36-T2	LL2-AEV42-T2	LL2-AEV48-T2	LL2-AEV54-T2
	GH : 2.0	-	LL2-AEV3620-T2	LL2-AEV4220-T2	LL2-AEV4820-T2	LL2-AEV5420-T2
LL-Multibase (Bridge)	GH : 1.3	LL2-AEV30-N2	LL2-AEV36-N2	LL2-AEV42-N2	LL2-AEV48-N2	LL2-AEV54-N2
	GH : 2.0	-	LL2-AEV3620-N2	LL2-AEV4220-N2	LL2-AEV4820-N2	LL2-AEV5420-N2
LS-Multibase / Single		-	-	-	-	-
CMFit Abutment / Ø10		CMFit10-AEV30-T2	CMFit10-AEV36-T2	CMFit10-AEV42-T2	CMFit10-AEV48-T2	CMFit10-AEV54-T2
CMFit Abutment / Ø14		-	CMFit14-AEV36-T2	CMFit14-AEV42-T2	CMFit14-AEV48-T2	CMFit14-AEV54-T2
3D-Analog		3D-AEV30	3D-AEV36	3D-AEV42	3D-AEV48	3D-AEV54
Screw(Sea)		AC-AEV305	AC-AEV365	AC-AEV425	AC-AEV48545	AC-AEV48545

Compatible Implant Systems

Camlog


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.28Hex
				20Ncm

Description-Library code	1401- CAM33	1402- CAM38	1403- CAM43	1404- CAM50	1405- CAM60
Connection Platform Ø	 Ø3.3	 Ø3.8	 Ø4.3	 Ø5.0	 Ø6.0
Product					
Scanbody	SC-CAM33	SC-CAM38	SC-CAM43	SC-CAM5060	SC-CAM5060
LL-Multibase / Single	LL2-CAM33-T2	LL2-CAM38-T2	LL2-CAM43-T2	LL2-CAM50-T2	LL2-CAM60-T2
LL-Multibase / Bridge	LL2-CAM33-N2	LL2-CAM38-N2	LL2-CAM43-N2	LL2-CAM50-N2	LL2-CAM60-N2
LS-Multibase / Single	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-CAM33-T2	CMFit10-CAM38-T2	CMFit10-CAM43-T2	CMFit10-CAM50-T2	CMFit10-CAM60-T2
CMFit Abutment / Ø14	-	CMFit14-CAM38-T2	CMFit14-CAM43-T2	CMFit14-CAM50-T2	CMFit14-CAM60-T2
3D-Analog	3D-CAM33	3D-CAM38	3D-CAM43	3D-CAM5060	3D-CAM5060
Screw(Sea)	AC-CAM5	AC-CAM5	AC-CAM5	AC-CAM50605	AC-CAM50605

Conelog

*Set (Scanbody/Multibase/CMFit Abutment + screw)

GMS-Scanbody (Intra-Oral)	Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
 (11.5mm)					 1.28Hex
					20Ncm

Description-Library code	1406- CN33	1407- CN3843	1408- CN50
Connection Platform Ø	 Ø3.3	 Ø3.8/4.3	 Ø5.0/6.0
Product			
GMS-Scanbody(Long)	GMS-CN33L	GMS-CN3843L	GMS-CN50L
CMFit Abutment / Ø10	CMFit10-CN33-E	-	-
CMFit Abutment / Ø14	-	CMFit14-CN3843-E	CMFit14-CN50-E
3D-Analog	3D-CN33	3D-CN3843	3D-CN50
Screw(Sea)	AC-CN5	AC-CN5	AC-CN505

Compatible Implant Systems

MIS SEVEN

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description-Library code	1501 - MS33	1502 - MSR	1503 - MSW
Connection Platform Ø	Narrow Ø3.3	Standard Ø 3.75 / 4.2	Wide Ø 5.0 / 6.0
Product			
SC-Scanbody	SC-MS33	SC-TSV3	SC-TSV4
GMS-Scanbody(Short)	-	GMS-TSV3S	GMS-TSV4S
GMS-Scanbody(Long)	-	GMS-TSV3L	GMS-TSV4L
LL-Multibase / Single	LL2-MS33-H	LL2-MSR-H	LL2-MSW-H
LL-Multibase / Bridge	LL2-MS33-N	LL2-MSR-N	LL2-MSW-N
LS-Multibase / Single	-	-	-
CMFit Abutment / Ø10	CMFit10-MS33-H	CMFit10-MSR-H	CMFit10-MSW-H
CMFit Abutment / Ø14	-	CMFit14-MSR-H	CMFit14-MSW-H
3D-Analog	3D-MS33	3D-TSV3	3D-TSV4
Screw(Sea)	AC-MS335	AC-MS5	AC-MS5

Compatible Implant Systems

MIS C1

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
				-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description—Library code	1504- C1M	1505- C1R	1506- C1W
Connection Platform Ø	C1  Ø3.3	C1  Ø 3.75/4.2	C1  Ø 5.0
Product			
Scanbody	SC-C1M	SC-C1R	SC-C1W
GMS-Scanbody(Short)	GMS-C1MS	GMS-C1RS	GMS-C1WS
GMS-Scanbody(Long)	GMS-C1ML	GMS-C1RL	GMS-C1WL
LL-Multibase / Single	LL2-C1M-T	LL2-C1R-T	LL2-C1W-T
LL-Multibase / Bridge	LL2-C1M-N	LL2-C1R-N	LL2-C1W-N
LS-Multibase / Single	-	-	-
CMFit Abutment / Ø10	CMFit10-C1M-T	CMFit10-C1R-T	CMFit10-C1W-T
CMFit Abutment / Ø14	CMFit14-C1M-T	CMFit14-C1R-T	CMFit14-C1W-T
3D-Analog	3D-C1M	3D-C1R	3D-C1W
Screw(Sea)	AC-C1M5	AC-C15	AC-C15

Compatible Implant Systems

SIC ace / max

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.2Hex
				20Ncm

SIC ace / max

Description-Library code	1601- SICM	1602- SICR
Connection Platform Ø	 Ø 3.3	 Ø 4.2
Product		
Scanbody	SC-SIC	SC-SIC
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	-	-
LL-Multibase / Single	LL2-SICM-H	LL2-SICR-H
LL-Multibase / Bridge	LL2-SICM-N	LL2-SICR-N
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-SICM-H	CMFit10-SICR-H
CMFit Abutment / Ø14	CMFit14-SICM-H	CMFit14-SICR-H
3D-Analog	3D-SICM	3D-SICR
Screw(5ea)	AC-SIC5	AC-SIC5

Compatible Implant Systems

Thommen SPI

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Four lobe 15Ncm Ø3.5 25Ncm Ø 4.0/4.5/5.0/6.0

Description—Library code	1701- SPI35	1702- SPI40	1703- SPI45	1704- SPI50	1705- SPI60
Product	Connection PlatformØ  Ø3.5	 Ø4.0	 Ø4.5	 Ø5.0	 Ø6.0
Scanbody	SC-SPI35	SC-SPI40	SC-SPI45	SC-SPI45	SC-SPI45
GMS-Scanbody(Short)	-	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-	-
LL-Multibase / Single	LL2-SPI35-H	LL2-SPI40-H	LL2-SPI45-H	LL2-SPI50-H	LL2-SPI60-H
LL-Multibase / Bridge	LL2-SPI35-N	LL2-SPI40-N	LL2-SPI45-N	LL2-SPI50-N	LL2-SPI60-N
LS-Multibase / Single	-	-	-	-	-
LS-Multibase / Bridge	-	-	-	-	-
CMFit Abutment / Ø10	CMFit10-SPI35-H	CMFit10-SPI40-H	CMFit10-SPI45-H	CMFit10-SPI50-H	CMFit10-SPI60-H
CMFit Abutment / Ø14	-	CMFit14-SPI40-H	CMFit14-SPI45-H	CMFit14-SPI50-H	CMFit14-SPI60-H
3D-Analog	3D-SPI35	3D-SPI40	3D-SPI45	3D-SPI50	-
Screw(Sea)	AC-SPI355	AC-SPI5	AC-SPI5	AC-SPI5	AC-SPI5

Compatible Implant Systems

BEGO

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25Hex
				30Ncm

Description-Library code	1801 - BEG037	1802 - BEG041	1803 - BEG045	1804 - BEG055
Connection Platform Ø	 Ø3.7	 Ø4.1	 Ø4.5	 Ø5.5
Product				
SC-Scanbody	SC-BEG0	SC-BEG0	SC-BEG0	SC-BEG0
GMS-Scanbody(Short)	-	-	-	-
GMS-Scanbody(Long)	-	-	-	-
LL-Multibase / Single	LL2-BEG037-H	LL2-BEG041-H	LL2-BEG045-H	LL2-BEG055-H
LL-Multibase / Bridge	LL2-BEG037-N	LL2-BEG041-N	-	-
LS-Multibase / Single	-	-	-	-
CMFit Abutment / Ø10	CMFit10-BEG037-H	CMFit10-BEG041-H	CMFit10-BEG045-H	CMFit10-BEG055-H
CMFit Abutment / Ø14	-	-	CMFit14-BEG045-H	CMFit14-BEG055-H
3D-Analog	3D-BEG037	3D-BEG041	3D-BEG045	3D-BEG055
Screw(Sea)	AC-BEG05	AC-BEG05	AC-BEG05	AC-BEG05

Compatible Implant Systems

MYTIS Arrow Implant

*set (Scanbody/Link/CMFit Abutment + screw)

OS-Scanbody(Intra Oral)	GMS-Scanbody(Intra-Oral)	Ø10 CMFit Abutment	Screw	Driver (Ncm)
	 (11.5mm)			 1.3Hex 30Ncm

Description-Library code	1901 - AR3740	1902 - AR4760
Connection Platform Ø	 Ø 3.7/4.0	 Ø 4.7/6.0
Product		
OS-Scanbody	-	OS-AR4760
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	GMS-AR3740L	-
CMFit Abutment / Ø10	CMFit10-AR3740-H	CMFit10-AR4760-H
Screw(5ea)	AC-AR5	AC-AR5

Compatible Implant Systems

Kyocera POI


*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	Single LL-Multibase	Bridge LL-Multibase	Ø10 CMFit Abutment	Screw	Driver (Ncm)
		-			 1.1Hex 20Ncm

Description-Library code	2001 - POI37	2002 - POI42
Connection Platform Ø Product	 Ø3.7	 Ø4.2
SC-Scanbody	SC-POI37	SC-POI42
LL-Multibase / Single	LL2-POI37-H	LL2-POI42-H
LL-Multibase / Bridge	-	-
CMFit Abutment / Ø10	CMFit10-POI37-H	CMFit10-POI42-H
Screw(5ea)	AC-POI375	AC-POI425

Kyocera FINESIA

*set (Scanbody/Link/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 Torx 20Ncm

Description-Library code	3601- FBLN	3602- FBLR
Connection Platform Ø Product	 Ø 3.2 / 3.4	 Ø 3.7 / 4.2
GMS-Scanbody(Short)	-	-
GMS-Scanbody(Long)	GMS-FBLNL	GMS-FBLRL
LL-Multibase / Single	-	LL2-FBLR-H
LL-Multibase / Bridge	-	-
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-FBLN-H	CMFit10-FBLR-H
CMFit Abutment / Ø14	-	-
3D-Analog	3D-FBLN	3D-FBLR
Screw(5ea)	AC-FBLN5	AC-FBLR5

Compatible Implant Systems

Biohorizons Laser Lok 3.0

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.25Hex
				30Ncm

Description-Library code	2101 - BH30
Connection Platform Ø	 Ø 3.0
Product	
OS-Scanbody	OS-BH30
LL-Multibase / Single	LL2-BH30-H2
LL-Multibase / Bridge	LL2-BH30-N2
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-BH30-H2
CMFit Abutment / Ø14	-
3D-Analog	3D-BH30
Screw(5ea)	AC-BH305

Compatible Implant Systems

Anthogyr AXIOM

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm) (7.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.18Hex 25Ncm

Description-Library code	2201 - ANG
Connection Platform Ø	universal  Ø 3.4/4.0/4.6/5.2
Product	
SC-Scanbody	SC-ANG
GMS-Scanbody(Short)	GMS-ANGS
GMS-Scanbody(Long)	GMS-ANGL
LL-Multibase / Single	LL2-ANG-E
LL-Multibase / Bridge	LL2-ANG-N
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-ANG-E
CMFit Abutment / Ø14	-
3D-Analog	3D-ANG
Screw(5ea)	AC-ANG5

Compatible Implant Systems

Warantec Oneplant


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Ø10 CMFit Abutment	Single LS-Multibase	Screw	Driver (Ncm)
		-		1.2Hex 30Ncm

Description-Library code	2301 - WOP		
Connection Platform Ø	Standard		Ø 4.3/4.8/5.3
Product			
Scanbody	SC-WOP		
LS-Multibase / Single	-		
CMFit Abutment Ø10	CMFit10-WOP-0		
Screw(5ea)	AC-WOP5		

Warantec IU

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			1.2Hex 15Ncm Ø 3.5/4.0/4.5/5.5

Description-Library code	2302 - SUMW		
Connection Platform Ø	Mini		Ø 3.3/3.6
Product			
SC-Scanbody	-		
GMS-Scanbody(Short)	-		
GMS-Scanbody(Long)	GMS-SUMWL		
LL-Multibase / Single	LL2-SUMW-H2		
CMFit Abutment / Ø10	CMFit10-SUMW-H2		
CMFit Abutment / Ø14	-		
3D-Analog	3D-SUMW		
Screw(5ea)	AC-SUMW5		

Compatible Implant Systems

Tekka In-Kone

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.2Hex
				15Ncm Ø 3.5/4.0/4.5/5.5

Description-Library code	2401 - TEK		
Connection Platform Ø	Standard		Ø 3.5/4.0/4.5/5.5
Product			
SC-Scanbody	SC-TEK		
GMS-Scanbody(Short)	-		
GMS-Scanbody(Long)	-		
LS-Multibase / Single	-		
CMFit Abutment / Ø10	CMFit10-TEK-H		
CMFit Abutment / Ø14	-		
3D-Analog	3D-TEK		
Screw(5ea)	AC-TEK5		

Compatible Implant Systems

Keystone PrimaConnex / Max-TL

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.25 Square 30Ncm

Description—Library code	2501 – KEP35	2502 – KEP41	2503 – KEP50	2504 – KEX50	2505 – KEX60
Connection Platform Ø	Ø 3.5	Ø 4.1	Ø 5.0 / 6.0	Ø 7.0	Ø 8.0
Product					
SC-Scanbody	SC-KEP35	SC-KEP41	SC-KEP50	SC-KEP50	SC-KEP50
GMS-Scanbody(Short)	GMS-KEP35S	GMS-KEP41S	GMS-KEP50S	GMS-KEP50S	GMS-KEP50S
GMS-Scanbody(Long)	GMS-KEP35L	GMS-KEP41L	GMS-KEP50L	GMS-KEP50L	GMS-KEP50L
LL-Multibase / Single	LL2-KEP35-T2	LL2-KEP41-T2	LL2-KEP50-T2	LL2-KEX50-T2	LL2-KEX65-T2
LL-Multibase / Bridge	LL2-KEP35-N2	LL2-KEP41-N2	LL2-KEP50-N2	LL2-KEX50-N2	LL2-KEX65-N2
LS-Multibase / Single	LS2-KEP35-T2	LS2-KEP41-T2	LS2-KEP50-T2	LS2-KEX50-T2	LS2-KEX65-T2
CMFit Abutment / Ø10	CMFit10-KEP35-T2	CMFit10-KEP41-T2	CMFit10-KEP50-T2	-	-
CMFit Abutment / Ø14	CMFit14-KEP35-T2	CMFit14-KEP41-T2	CMFit14-KEP50-T2	CMFit14-KEX50-T2	CMFit14-KEX65-T2
3D-Analog	3D-KEP35	3D-KEP41	3D-KEP50	3D-KEX50	3D-KEX65
Screw(Sea)	AC-KEPT5	AC-KEPT5	AC-KEPT5	AC-KEPT5	AC-KEPT5

Compatible Implant Systems

Neodent CM


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 0.9Hex 15Ncm

Description-Library code	2601 - NCM			
Product	Connection Platform Ø	Standard		Ø 3.5/4.3/5.0
SC-Scanbody	SC-NCM			
LS-Multibase / Single	-			
CMFit Abutment / Ø10	CMFit10-NCM-H			
CMFit Abutment / Ø14	-			
3D-Analog	3D-NCM			
Screw	-			

Neodent GM

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	 (11.5mm)			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 unigrip (not compatible with Nobel) 20Ncm

Description—Library code		2602 - NCG		
Product	Connection Platform Ø	Standard		Ø 3.5/3.75/4.0/4.3/5.0
SC—Scanbody		SC-NCG		
GMS—Scanbody(Long)		GMS-NCGL		
CMFit Abutment / Ø10		CMFit10-NCG-H2		
CMFit Abutment / Ø14		CMFit14-NCG-H2		
LL—Multibase (Single)	GH : 1.3	LL2-NCG-H2		
	GH : 2.0	LL2-NCG20-H2		
LL—Multibase (Bridge)	GH : 1.3	LL2-NCG-N2		
	GH : 2.0	LL2-NCG20-N2		
3D—Analog		3D-NCG		
Screw(5ea)		AC-NCG5		

Neodent CM

Neodent GM

Compatible Implant Systems

Bicon Implant

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Bridge LL-Multibase	Bridge LS-Multibase	Ø10 CMFit Abutment	Ø14 CMFit Abutment
		-		-

Description-Library code	2701- BC20	2702- BC25	2703- BC30
Connection Platform Ø	○ Ø 2.0	○ Ø 2.5	○ Ø 3.0
Product			
SC-Scanbody	SC-BC20	SC-BC25	SC-BC30
LL-Multibase / Bridge	LL2-BC20-N	LL2-BC25-N	LL2-BC30-N
LS-Multibase / Bridge	-	-	-
CMFit Abutment / Ø10	CMFit10-BC20-N	CMFit10-BC25-N	CMFit10-BC30-N
CMFit Abutment / Ø14	-	-	-

Compatible Implant Systems

Medentika IPS

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Ø10 CMFit Abutment	Ø14 CMFit Abutment	Screw	Driver (Ncm)
		-		 1.3 Hex 25Ncm

Description-Library code	2801- MTR
Connection Platform Ø	 Ø3.5/4.3/5.0
Product	
SC-Scanbody	SC-MTR
CMFit Abutment / Ø10	CMFit10-MTR-R
Screw(5ea)	AC-MTR5

GC Aadva Implant

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	Ø10 CMFit Abutment	Ø14 CMFit Abutment	Screw	Driver (Ncm)
		-		 Torx 20Ncm

Description-Library code	2901- GAD
Connection Platform Ø	Regular / Wide  Ø4.0 / 5.0
Product	
SC-Scanbody	SC-GAD
CMFit Abutment / Ø10	CMFit10-GAD-H
Screw(5ea)	AC-GAD5

Medentika IPS

GC Aadva Implant

Compatible Implant Systems

BIOTECH KONTACT

*Set (Scanbody/Multibase/CMFit Abutment + screw)



SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.2 Hex 20Ncm

Description-Library code	3001- BK30	3002- BK40
Connection Platform Ø  Ø 3.0		 Ø 3.6 / 4.2 / 4.8 / 5.4
Product		
SC-Scanbody	SC-BK30	SC-BK40
LL-Multibase / Single	LL2-BK30-T	LL2-BK40-T
LL-Multibase / Bridge	LL2-BK30-N	LL2-BK40-N
LS-Multibase / Single	-	-
CMFit Abutment / Ø10	CMFit10-BK30-T	CMFit10-BK40-T
3D-Analog	3D-BK30	3D-BK40
Screw(Sea)	AC-BK305	AC-BK405

Compatible Implant Systems

Medentis medical ICX-System

*Set (Scanbody/Multibase/CMFit Abutment + screw)


SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
	-			-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 1.42Hex 30Ncm

Description-Library code	3101 - ICX
Connection Platform Ø	universal  Ø3.75/4.1/4.8
Product	
SC-Scanbody	SC-ICX
LL-Multibase / Single	LL2-ICX-H
LL-Multibase / Bridge	LL2-ICX-N
LS-Multibase / Single	-
CMFit Abutment / Ø10	CMFit10-ICX-H
CMFit Abutment / Ø14	CMFit14-ICX-H
3D-Analog	3D-ICX
Screw(5ea)	AC-ICX5

Compatible Implant Systems

B&B Dental


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
	-			 1.27 Hex 25Ncm

Description—Library code	3201- BNB
Product	Connection Platform Ø  Ø 4.0/4.5/5.0
GMS-Scanbody(long)	GMS-BNBL
LL-Multibase / Single	LL2-BNB-H
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	CMFit10-BNB-H
3D-Analog	3D-BNB
Screw(5ea)	AC-BNB5

Bredent Medical Sky

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
				 Torx 25Ncm

Description—Library code	3301- SKR
Product	Connection Platform Ø  Ø 3.5(NP) / 4.0-6.0
GMS-Scanbody(long)	GMS-SKRL
LL-Multibase / Single	LL2-SKR-T
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	CMFit10-SKR-T
CMFit Abutment / Ø14	CMFit14-SKR-T
3D-Analog	3D-SKR
Screw(5ea)	AC-SKR5

Compatible Implant Systems

NEOSS ProActive


*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)		-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-				 Unigrip (not compatible with Nobel)
				32Ncm

Description—Library code	3401- NEP
Connection Platform Ø	 Ø3.5/4.0/4.5/5.0/5.5
Product	
GMS-Scanbody(long)	GMS-NEPL
LL-Multibase / Single	LL2-NEP-T
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	-
CMFit Abutment / Ø14	CMFit14-NEP-T
3D-Analog	3D-NEP
Screw(5ea)	AC-NEP5

Champions (R)evolution

*Set (Scanbody/Multibase/CMFit Abutment + screw)

SC-Scanbody(Model)	GMS-Scanbody(Intra-Oral)	Single LL-Multibase	Bridge LL-Multibase	Single LS-Multibase
-	 (11.5mm)	-	-	-
Ø10 CMFit Abutment	Ø14 CMFit Abutment	3D-Analog	Screw	Driver (Ncm)
-				 1.2 Hex
				30Ncm

Description—Library code	3501- CHR
Connection Platform Ø	 Ø 3.5/4.5/5.5
Product	
GMS-Scanbody(long)	GMS-CHRL
LL-Multibase / Single	-
LL-Multibase / Bridge	-
CMFit Abutment / Ø10	-
CMFit Abutment / Ø14	CMFit14-CHR-H
3D-Analog	3D-CHR
Screw(5ea)	AC-CHR5

Compatible Implant Systems

MultiUnit(Multibase)

Dentium Screw Abutment

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
		-		 1.27Hex
				10Ncm

Description-Library code	0203- MUDE45	0204- MUDE55
Connection Platform Ø	Ø 4.5	Ø 5.5
Product		
SC-Scanbody	SC-MUDE	SC-MUDE
LB-Multibase / Bridge	LB2-MUDE45	LB2-MUDE55
3D Analog	-	-
Screw(Sea)	AC-MUDE5	AC-MUDE5

Straumann Tissue Level SynOcta coping

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
		-		 Torx
				35Ncm

Description-Library code	0706- MUIT
Connection Platform Ø	Ø 4.8
Product	
SC-Scanbody	SC-MUIT
LB-Multibase / Bridge	LB2-MUIT
3D Analog	-
Screw(Sea)	AC-MUIT5

Straumann Bonelevel SRA

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 Torx
				15Ncm

Description-Library code	0707- MUSR35	0708- MUSR46
Connection Platform Ø	Ø 3.5	Ø 4.6
Product		
SC-Scanbody	SC-MUSR35	SC-MUSR46
LB-Multibase / Bridge	LB2-MUSR35	LB2-MUSR46
3D Analog	3D-MUSR35	3D-MUSR46
Screw(Sea)	AC-MUSR5	AC-MUSR5

Compatible Implant Systems

MultiUnit (Multibase)

Nobel MultiUnit 4.8

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 unigrip
				15Ncm

Description–Library code	0812- MUTB
Product	Connection Platform Ø Ø 4.8
SC-Scanbody	SC-MUTB
LB-Multibase / Bridge(4.0mm)	LB2-MUTB2
LB-Multibase / Bridge(5.5mm)	LB2-MUTBL2
3D Analog	3D-MUTB
Screw(5ea)	AC-MUTB5

Ankylos Multi-Unit

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Link	3D-Analog	Screw	Driver (Ncm)
		-		 1.0Hex
				10Ncm

Description–Library code	0907- AKLM	0908- AKLR
Product	Connection Platform Ø Mini	Regular
SC-Scanbody	SCB-16	SCB-16
LB-Multibase / Bridge	LB2-AKLM	LB2-AKLR
3D Analog	-	-
Screw(5ea)	AC-AKL5	AC-AKL5

Zimmer Tapered Screw-vent Multi-Unit(TSV)

*set (Scanbody/Multibase + screw)

SC-Scanbody	Bridge LB-Multibase	3D-Analog	Screw	Driver (Ncm)
				 1.27Hex
				20Ncm


Description–Library code	1206- MUTSV
Product	Connection Platform Ø Universel
SC-Scanbody	SC-MUTSV
LB-Multibase / Bridge	LB2-MUTSV
3D Analog	3D-MUTSV
Screw(5ea)	AC-MUTSV5

GTS Transfer Scanbody

GTS Transfer Scanbody

Digital implant technology

Convenient & prompt digital implant impression
Ideal subgingiva designing referenced from digital library



Differentiate the Gingival height by appearance



- Biocompatible **Titanium (Ti-6Al-4V ELI)** from **CARPENTER Co.** guarantees high-durability.
- Autoclavable.
- **Providing various sizes.**
- Osstem TS And Dentium Super Line Connections are available.
- Providing **High precise Geo library** for CAD s/w(3Shape, Exocad, Dental wings)
- **Special tool** allows finding right position and preventing detachment of screw.
- **Zirconia prosthesis designing** with multi-base

Guide Tool



Guide Tool protects the scan body from any damage and helps to find dental implant connection easily.

Healing Abutment & Bite Index



GTS Scanbody is multifunctional that can be also used as a healing abutment and impression coping.


The groove on the scan body allows users to define different gingiva types and provide a bite index when taking an impression.

GTS Transfer scanbody


Direct Impression → Digital process



GTS scanbody in patient's mouth



Impression taking and stone model casting



Designing process with scan data from model scanner



Designing process with scan data from model scanner



Final prosthesis setting to the patient's mouth

GTS Transfer Scanbody List



Guide Tool



DER – Series (Dentium Super-line / Megagen AnyOne / Dentis S-clean / Neo IS)

Diameter GH / Height	4.0mm	4.5mm	5.5mm	6.5mm
 2.0mm / 5.0mm	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
 4.0mm / 7.0mm	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
 6.0mm / 9.0mm	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560
Guide Tool	GTST-D40-127H	GTST-D45-127H	GTST-D55-127H	GTST-D65-127H

* Laser mark [DER55: Diameter 5.5mm]
Easy verifying of gingival heights by differentiated appearance.

SUM – Series (Osstem Implant TS Mini)

Diameter GH / Height	4.0mm	4.5mm
 2.0mm / 5.0mm	GTS-SUM4020	GTS-SUM4520
 4.0mm / 7.0mm	GTS-SUM4040	GTS-SUM4540
 6.0mm / 9.0mm	GTS-SUM4060	GTS-SUM4560
Guide Tool	GTST-S40-120H	GTST-S45-120H

* Laser mark [SUM45: Diameter 4.5mm]
Easy verifying of gingival heights by differentiated appearance.

SURO – Series (Osstem Implant TS Standard / Dio UF Regular)

Diameter GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
 2.0mm / 5.0mm	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
 4.0mm / 7.0mm	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
 6.0mm / 9.0mm	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
Guide Tool	GTST-S45-120H	GTST-S50-120H	GTST-S60-120H	GTST-S70-120H

* Laser mark [SURO50: Diameter 5.0mm]
Easy verifying of gingival heights by differentiated appearance.

Sirona Ti-base

Sirona Multibase



Scan Cap for Sirona Multibase

Scan Cap

Scan Cap (short)	Scan Cap (long)
SCSN-CAP	SCSNL-CAP

Description

Osstem Implant TS

*Set (Sirona Multibase + screw)

 Single	Platform Ø	TS Mini	TS Standard	Driver (Ncm)	1.2Hex
	Connection	Ø3.0/3.5	Ø4.0/4.5/5.0		
	Single [GH:1.2mm]	SRN2-SUM-H2	SRN2-SUR0-H2	20Ncm	Ø3.0/3.5
	Single [GH:2.6mm]	-	SRN2-SUR026-H2	30Ncm	Ø4.0/4.5/5.0

Osstem Implant US / Biomet 3i Osseotite 4.1

*Set (Sirona Multibase + screw)

 Single	Platform Ø	US Regular	Driver (Ncm)	1.2Hex
	Connection	Ø4.1		
	Single	SRN2-EXR-H2	30Ncm	

Sirona Ti-base

Sirona Multibase

Description

Dentium Implantium / SuperLine

*Set (Sirona Multibase + screw)

 Single	Platform Ø	universal  Ø3.7/4.0/4.5		Driver (Ncm)  1.27Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DER-H2		25-30Ncm
Single [GH:2.6mm]	SRN2-DER26-H2			

Megagen AnyRidge

*Set (Sirona Multibase + screw)

 Single	Platform Ø	Universal  Ø3.5/4.0/4.5		Driver (Ncm)  1.2Hex
	Connection			
	Single	SRN2-MARR-H2		35Ncm

Megagen AnyOne

*Set (Sirona Multibase + screw)

 Single	Platform Ø	Universal  Ø 3.5/4.0/4.5/5.0		Driver (Ncm)  1.2Hex
	Connection			
	Single [GH:1.2mm]	SRN2-DERA-H2		25Ncm
Single [GH:2.6mm]	SRN2-DERA26-H2			

Neo Biotech IS

*Set (Sirona Multibase + screw)

 Single	Platform Ø	 Ø3.5/4.0/4.5	Driver (Ncm)  1.2Hex
	Connection		
	Single [GH:1.2mm]	SRN2-DERN-H2	30Ncm
Single [GH:2.6mm]	SRN2-DERN26-H2		

Straumann Bone Level

*Set (Sirona Multibase + screw)

 Single	Platform Ø	NC  Ø3.3	RC  Ø4.1/4.8	Driver (Ncm)  Torx
	Connection			
	Single	SRN2-BLM-R2	SRN2-BLR-R2	35Ncm

Sirona Ti-base

Sirona Multibase

Description

Straumann Tissue Level

*Set (Sirona Multibase + screw)

 Single	Platform Ø	RN  Ø4.8	WN  Ø6.5	Driver (Ncm)  Torx
	Connection	SRN2-ITR-02	SRN2-ITW-02	35Ncm

Nobel Biocare Active

*Set (Sirona Multibase + screw)

 Single	Platform Ø	NP  Ø3.5	RP  Ø4.3/5.0	Driver (Ncm)  unigrip
	Connection	SRN2-ACN-H2	SRN2-ACR-H2	35Ncm

Nobel Biocare Branemark

*Set (Sirona Multibase + screw)

 Single	Platform Ø	RP  Ø4.1	Driver (Ncm)  unigrip
	Connection	SRN2-EXRB-H2	35Ncm

Nobel Biocare Replace

*Set (Sirona Multibase + screw)

 Single	Platform Ø	NP  Ø3.5	RP  Ø4.3	WP  Ø5.0	Driver (Ncm)  unigrip
	Connection	SRN2-RPM-E2	SRN2-RPR-E2	SRN2-RPW-E2	35Ncm

Dentsply Xive / Frialit2

*Set (Sirona Multibase + screw)

 Single	Platform Ø	 Ø3.4	 Ø3.8	 Ø4.5	Driver (Ncm)  1.22Hex
	Connection	SRN2-XVM-H2	SRN2-XVR-H2	SRN2-XVW-H2	24Ncm

Sirona Ti-base

Sirona Multibase

Description


Dentsply Ankylos

*Set (Sirona Multibase + screw)

 Single Bridge	Platform Ø	universal  Ø 3.4/4.5/5.5		Driver (Ncm)	15Ncm
	Connection			 1.0Hex	
	Single [GH:1.0mm]	SRN2-AKLL1-T			
	Single [GH:2.0mm]	SRN2-AKLL2-T			
	Bridge [GH:1.0mm]	SRN2-AKLL1-N			
	Bridge [GH:2.0mm]	SRN2-AKLL2-N			

Zimmer Tapered Screw-vent

*Set (Sirona Multibase + screw)

 Single	Platform Ø	 Ø3.5	 Ø4.5	 Ø5.7	Driver (Ncm)
	Connection				 1.25Hex
	Single	SRN2-TSV3-H2	SRN2-TSV4-H2	SRN2-TSV5-H2	30Ncm

Astra Tech Osseospeed

*Set (Sirona Multibase + screw)

 Single	Platform Ø	Aqua  Ø 3.5/4.0		Lilac  Ø 4.5/5.0		Driver (Ncm)
	Connection					 1.27Hex
	Single	SRN2-ASR-H2		SRN2-ASW-H2		20Ncm Ø 3.5/4.0
						25Ncm Ø 4.5/5.0

CCM Casting Abutment

Product Line

CCM Casting Abutment

Reliable product with accumulated clinical data since 2009
Excellent bio-compatibility
Optimized design for emergence profile.
All type of implant restorations are enable.



Description

Osstem Implant TS

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø Connection	TS Mini  Ø3.0/3.5	TS Standard  Ø4.0/4.5/5.0		Driver (Ncm)  1.2Hex	
Single	SUM-H2		GH : 1mm	SUR0-H2	20Ncm	Ø3.0/3.5
			GH : 2.2mm	SUR022-H2		
			GH : 3mm	SUR030-H2		
Bridge	SUM-N2		GH : 1mm	SUR0-N2	30Ncm	Ø4.0/4.5/5.0
			GH : 2.2mm	SUR022-N2		
			GH : 3mm	SUR030-N2		

Osstem Implant US

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø Connection	US Mini  Ø3.5	US Regular  Ø4.1	US Wide  Ø5.1	Driver (Ncm)  1.2Hex	
Single	EXM-H2	EXR-H2	EXW-H2		30Ncm	
Bridge	EXM-N2	EXR-N2	EXW-N2			

Osstem Implant SS

*Set (CCM Abt + screw)


 Single Bridge	Platform Ø Connection	SS Regular  Ø4.8	SS Wide  Ø6.0	Driver (Ncm)  1.2Hex	
Single	ISR-02	ISW-02		30Ncm	
Bridge	ISR-N2	ISW-N2			

CCM Casting Abutment

Description


Dentium SuperLine / Implantium

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	Standard  Ø3.7/4.0/4.5		Driver (Ncm)  1.27Hex
	Single	GH : 1mm	DER-H2	25-30Ncm
		GH : 2mm	DER20-H2	
		GH : 3mm	DER30-H2	
	Bridge	GH : 1mm	DER-N2	
		GH : 2mm	DER20-N2	
		GH : 3mm	DER30-N2	


Megagen Anyridge

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	Standard  Ø3.5/4.0/4.5	Driver (Ncm)  1.2Hex
	Single	MARR-H2	35Ncm
	Bridge	MARR-N2	

Megagen Ex-feel [Internal]

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø 4.8	Driver (Ncm)  1.2Hex
	Single	ITRK-02	35Ncm
	Bridge	ITRK-N2	



- Qualified Raw Material
Cobalt-chrome (CoCr) alloys
- Single and bridge case
Available single and bridge types of abutments
- Two screws
TiN coating screw & Titanium screw
- Compatibility
Various compatible Implant systems

CCM Casting Abutment

Description

DIO UF / UFII

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø 3.8/4.0/4.5/5.0	Driver (Ncm)  1.2Hex
	Connection		
	Single	SURD-H2	30Ncm
Bridge	SURD-N2		

Neo Biotech IS

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø		Ø3.5/4.0/4.5	Driver (Ncm)	 1.2Hex	
	Connection				30Ncm	
	Single	DERN-H2				
	Bridge	DERN-N2				

Straumann Bone Level

*Set (CCM Abt + screw)

 Single  Bridge  	Platform Ø Connection	NC  Ø3.3	RC  Ø4.1/4.8	Driver (Ncm)  Torx
	Single	BLM-R2	BLR-R2	35Ncm
	Bridge	BLM-N2	BLR-N2	

Straumann Tissue Level

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	RN	 Ø 4.8	WN	 Ø 6.5	Driver (Ncm)	 Torx
	Connection					35Ncm	
	Single	ITR-02		ITW-02			
Bridge	ITR-N2		ITW-N2				

Nobel Biocare Active

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	NP	 Ø3.5	RP	 Ø 4.3/5.0	Driver (Ncm)  unigrip
	Connection					35Ncm
	Single	ACN-H2		ACR-H2		
	Bridge	ACN-N2		ACR-N2		

Nobel Biocare Branemark

*Set (CCM Abt + screw)


    Single Bridge	Platform Ø Connection	NP  Ø3.5	RP  Ø4.1	WP  Ø5.1	Driver (Ncm)  unigrip
	Single	EXMB-H2	EXRB-H2	EXWB-H2	35Ncm
	Bridge	EXMB-N2	EXRB-N2	EXWB-N2	

CCM Casting Abutment

Description


Nobel Biocare Replace

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	NP  Ø3.5	RP  Ø4.3	WP  Ø5.0	Driver (Ncm)  unigrip
	Connection				
	Single	RPM-E2	RPR-E2	RPW-E2	35Ncm
	Bridge	RPM-N2	RPR-N2	RPW-N2	

Nobel Biocare Multi Unit

*Set (CCM Abt + screw)

 Bridge	Platform Ø	 Ø4.8			Driver (Ncm)  unigrip
	Connection				
	Bridge	MUTB-N2			15Ncm

Dentsply Xive / Frialit2

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø3.4	 Ø3.8	 Ø4.5	Driver (Ncm)  1.22Hex
	Connection				
	Single	XVM-H2	XVR-H2	XVW-H2	24Ncm
	Bridge	XVM-N2	XVR-N2	XVW-N2	


Dentsply Ankylos

*Set (CCM Abt + screw)

 Bridge	Platform Ø	A/B/C/D  Ø3.4/4.5/5.5	Driver (Ncm)  1.0Hex
	Connection		
	Single	-	25Ncm
Bridge	DNF3-N		

Biomet 3i Certain

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø3.4	 Ø4.1	Driver (Ncm)  1.2Hex
	Connection			
	Single	3ICM-H2	3ICR-H2	20Ncm
	Bridge	3ICM-N2	3ICR-N2	


Biomet 3i Osseotite

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø4.1	 Ø5.0	Driver (Ncm)  1.2Hex
	Connection			
	Single	EXR-H2	EXT-H2	20Ncm
	Bridge	EXR-N2	EXT-N2	

Zimmer Tapered Screw-vent

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	 Ø3.5	 Ø4.5	 Ø5.7	Driver (Ncm)  1.25Hex
	Connection				
	Single	TSV3-H2	TSV4-H2	TSV5-H2	30Ncm
	Bridge	TSV3-N2	TSV4-N2	TSV5-N2	

CCM Casting Abutment

Description


Zimmer Spline

*Set (CCM Abt + screw)

 Single	Platform Ø	 Ø 3.7 / 4.0		Driver (Ncm)	1.25Hex
	Connection	SP40-H		30Ncm	

Astra Tech Osseospeed

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	X-Small  Ø3.0	Aqua  Ø3.5/4.0	Lilac  Ø4.5/5.0	Driver (Ncm)	1.27Hex
	Connection	AS3-T2	ASR-H2	ASW-H2	15Ncm	Ø3.0
	Single	AS3-T2	ASR-H2	ASW-H2	20Ncm	Ø3.5/4.0
	Bridge	-	ASR-N2	ASW-N2	25Ncm	Ø4.5/5.0

MIS C1

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	C1  Ø3.3	C1  Ø3.75/4.2	C1  Ø5.0	Driver (Ncm)	1.25Hex
	Connection	C1M-T	C1R-T	C1W-T	30Ncm	
	Single	C1M-T	C1R-T	C1W-T		
	Bridge	C1M-N	C1R-N	C1W-N		

SIC ace / max

*Set (CCM Abt + screw)

 Single	Platform Ø	 Ø3.3	 Ø4.2	Driver (Ncm)	1.2Hex
	Connection	SICM-H	SICR-H	20Ncm	

Anthogyr AXIOM

*Set (CCM Abt + screw)

 Single Bridge	Platform Ø	universal  Ø 3.4/4.0/4.6/5.2	Driver (Ncm)	1.18Hex
	Connection			
	Single	ANG-E	25Ncm	
Bridge	ANG-N			

Tekka In-Kone

*Set (CCM Abt + screw)

 Single	Platform Ø	Standard  Ø3.5/4.0/4.5/5.5		Driver (Ncm)	1.2Hex
	Connection	TEK-H		15Ncm	

Geo Screw Driver



Angled Driver

Square Type / Ratchet Type

Use with LS Multibase Abutment

Maximum Angle 25°

Recommended torque allowance : below 30Ncm.(reference manufacture's guide)

Caution: unstable connection may bend or destroy the driver

Recommended usage (1.2hex/1.27hex)-10 times

(Universal Type)- 20 times

Connect with Square type Torque wrench

May use Round Type Torque wrench

by connecting HDEC - Short (Clinic)Handle



Connect Straight & Angle driver
with the Torque Wrench



Straight Driver

Square Type / Ratchet Type

Connect by manufacture's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench

by connecting HDEC - Short (Clinic)Handle

Universal - Round hexagon

Square Type / Ratchet Type

Connect by manufacture's recommended torque allowance

Damage on the tip or usage of long period of time need to exchange driver

Connect with Square type Torque wrench

May use Round Type Torque wrench

by connecting HDEC - Short (Clinic)Handle

Angled Driver

3i & Osstem - 1.2mm Hex

			
ASD3IL	ASD3IR	ALD3IL	ALD3IR
Square -20mm	Square -13mm	Ratchet -20mm	Ratchet -13mm

Astra & Zimmer & Dentium - 1.27mm Hex

			
ASDZIL	ASDZIR	ALDZIL	ALDZIR
Square -20mm	Square -13mm	Ratchet -20mm	Ratchet -13mm

Universal - Round hexagon

			
ASDUVL	ASDUVR	ALDUVL	ALDUVR
Square -20mm	Square -13mm	Ratchet -20mm	Ratchet -13mm

Geo Screw Driver

Straight Driver

3i - 1.2mm Hex

			
SSD3IL Square -20mm	SSD3IR Square -13mm	SLD3IL Ratchet -20mm	SLD3IR Ratchet -13mm
			
SSD3IEL Square -35mm			

Zimmer - 1.27mm Hex

			
SSDZIL Square -20mm	SSDZIR Square -13mm	SLDZIL Ratchet -20mm	SLDZIR Ratchet -13mm

Astra - 1.27mm Hex

			
SSDATL Square -20mm	SSDATR Square -13mm	SLDATL Ratchet -20mm	SLDATR Ratchet -13mm
			
SSDATEL Square -35mm			

Nobel - unigrip

			
SSDNBL Square -20mm	SSDNBR Square -13mm	SLDNBL Ratchet -20mm	SLDNBR Ratchet -13mm
			
SSDNBEL Square -35mm			

Straumann - Torx

			
SSDSTL Square -20mm	SSDSTR Square -13mm	SLDSTL Ratchet -20mm	SLDSTR Ratchet -13mm
			
SSDSTEL Square -35mm			

Geo Screw Driver

Driver Tool

 HDEL Long (Lab) Handle	 Long (Lab) Handle & Driver
 HDEC Short (Clinic) Handle	 Short (Clinic) Handle & Driver
 GTW Torque Wrench	 Round Type Torque Wrench & Handle Short & Driver  Square Type Torque Wrench & Driver

Lab Driver Set

		Description		
		Product		
		Size	Code	Implant System
		 1.2 Hex	Blue	Osstem Implant / Dio Implant Biomet 3i / Xive
		 1.27 Hex	Red	Dentium / Astra Tech Zimmer
		 Unigrip	Yellow	Nobel Biocare
		 Torx	Green	Straumann
		 1.0 Hex	Purple	Ankylos
		 Universal	Black	Dynamic



Lab Driver Set

Instructions of CAD/CAM Abutment

< Procedure for intra-oral scan >



1. Scanbody inserting



2. Scanning



3. Matching [Scanbody and Library]



4. Designing - 3shape / Exocad



5. Sintering



6. Setting

< Procedure for model scan >




1. Scanbody inserting






2. Scanning



3. Matching
[Scanbody and Library]



4. Designing - 3shape / Exocad



5. Sintering and model setting

Instructions of CCM Casting Abutment

< Screw Type manufacturing>



CCM Casting abutment



Status of Wax-up the connection of the metal part has more wax to prevent the occurrence of air bubble in the connection part.



Status after casting of general Porcelain metal.



Investment material and the oxide film is removed by using a Glass bead (50um/ 2~4 bar)



The oxide film can be removed cleanly without damaging the internal side of Hex.




The Junction is checked in the shape of work in the condition where the oxide film is removed.




Upon sanding with alumina oxide, it can be the cause of plaque.




Status after the oxide film is removed.



The formation of oxide film can be observed on the basal area of firing porcelain.



The firing porcelain part is wrapped by using wax.



The oxide film is removed by using the Glass bead.



The oxide film of internal side of Access hole is removed.



Final compatibility of finished upper structure is checked.




< Customized Type manufacturing>




Full contour wax up



Gum sculpture



after Casting



Instructions of CCM Casting Abutment

< Hybrid Type manufacturing>




Wax up on the Model



Separation of Abutment from Model



Sprue installation



After Casting



Screw Type is manufactured by using Denture Resin

< Bar Overdenture manufacturing>



Wax up on the Model






Milling in the desired form on the Wax



Bar Milling after Casting



Check the compliance



Free Geo Library available for 3Shape and exocad
High quality of stability and accuracy
One scanbody can be used to design all type of implant restorations
Providing library for Ti-base
Premill blank and implant bridge



Benefits :

Optimized to show best function in dental CAD
Enable to design individual prosthetic solutions for various implant systems.
Available in 3Shape and Exocad with various Ti-bases.
Flexible to fabricate milled structures with its individuality.



Application :

Individual ceramic abutments on titanium bases
Individual one-piece titanium abutments
Occlusal screwable crowns and bridges
Bars, telescopes
Screw retained abutment
Customized Zirconia abutment
Customized titanium abutment



Selection of Geo Library

Library for Implant System

Geo_0101_SUM
Geo_0102_SURO
Geo_0103_EXM
Geo_0104_EXR



Geo product code for different type

LS(short)Link
LL(long)Link
LL_C13_H40
LL_C13_H55
LL_C20_H43
LL_C20_H60
CMFit
Premill,blank

Gingival height(1.3mm/2.0mm)

Abutment height(cutting line)

Premill,blank



Cement gap option with correspond Scanbody

Cement [gap]-Zirconia

027 SC
035 SC
042 SC
052 SC
027 GMS-S
035 GMS-S
045 GMS-S
052 GMS-S
027 GMS-L
035 GMS-L
042 GMS-L
052 GMS-L
08 PMMA-SC
10 PMMA-SC
08 PMMA-GMS-S
10 PMMA-GMS-S
08 PMMA-GMS-L
10 PMMA-GMS-L



- Model Scanbody

- Intra-oral Scanbody

Height(7.0 ~8.0mm)

Height(10.5 ~11.5mm)

Cement [gap]-PMMA,Metal



Geo BS Library– Bridge&Single



Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
Osstem	TS M	SUM	✓	✓	SC-SUM	SUMC	v
	TS Standard	SURO	✓	✓	SC-SURO	SUROC	v
	US M	EXM	✓	✓	SC-EXM	EXMC	v
	US R	EXR	✓	✓	SC-EXR	EXRC	v
	US W	EXW	✓	✓	SC-EXW	EXWC	v
Dentium	Superline	DER	✓	✓	SC-DER	DERC	v
	Screw-Abutment	MUDE45	✓		SC-MUDE	MUDEC	v
		MUDE55	✓		SC-MUDE	MUDEC	v
Straumann	Bonelevel 3.3	BLM	✓	✓	SC-BLM	BLMC	v
	Bonelevel 4.1	BLR	✓	✓	SC-BLR	BLMC	v
	Bonelevel 4.8	BLR	✓	✓	SC-BLR	BLRC	v
	Tissue level 4.8	ITR	✓	✓	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	v
	Tissue level 6.5	ITW	✓	✓	SC-ITRW	bridge - ITTC	Temporary screw
						single - ITC	v
	SynOcta coping	MUIT	✓		SC-MUIT	MUITC	v
	SRA Abutment	MUSR35	✓		SC-MUSR35	MUSRC	v
		MUSR46	✓		SC-MUSR46	MUSRC	v
Nobel	Multi Unit	MUTB	✓		SC-MUTB	MUTBC	v
	Active 3.5	ACN	✓	✓	SC-ACN	ACNC	v
	Active 4.3	ACR	✓	✓	SC-ACR	ACRC	v
	Active 5.0	ACR	✓	✓	SC-ACR	ACRC	v
	Active 5.5	ACW	✓	✓	SC-ACW	ACRC	v
	Branemark 3.4	EXM	✓	✓	SC-EXM(B)	EXMBC	v
	Branemark 4.1	EXR	✓	✓	SC-EXR(B)	EXRBC	v
	Branemark 5.1	EXW	✓	✓	SC-EXW(B)	EXWBC	v
	Replace 3.5 rpm	RPM	✓	✓	SC-RPM	RPMC	v
	Replace 4.3 rpm	RPR	✓	✓	SC-RPR	RPRC	v
	Replace 5.0 rpm	RPW	✓	✓	SC-RPW	RPRC	v
	Replace 6.0 rpm	RP60	✓	✓	SC-RP60	RPRC	v
Dentsply	Xive 3.0	XV3	✓	✓	SC-XV3	XV3C	v
	Xive 3.4	XVM	✓	✓	SC-XVM	XVC	v
	Xive 3.8	XVR	✓	✓	SC-XVR	XVC	v
	Xive 4.5	XVW	✓	✓	SC-XVW	XVC	v
	Xive 5.5	XV55	✓	✓	SC-XVW	XVC	v
Ankylos	AK II	AKLL	✓		OS-AKLL	DNF3C	special screw

Geo BS Library- Bridge&Single

Brand	Implant System	System code	Bridge(Bar)	Single	Scanbody	Screw	
					code	Geo screw code	Compatible with original
3i	Certain 3.4	3ICM	✓	✓	SC-3ICM	3ICC	v
	Certain 4.1	3ICR	✓	✓	SC-3ICR	3ICC	v
	Certain 5.0	3ICW	✓	✓	SC-3ICR	3ICC	v
	Certain 6.0	3IC60	✓	✓	SC-3ICR	3ICC	v
	External 3.4	EX3M	✓	✓	SC-EX3M	EXRC	v
	External 4.1	EXR	✓	✓	SC-EXR	EXRC	v
	External 5.0	EXT	✓	✓	SC-EXR	EXRC	v
	External 6.0	EXTW	✓	✓	SC-EXR	EXRC	v
Zimmer	TSV 3.5	TSV3	✓	✓	SC-TSV3	TSVC	v
	TSV 4.5	TSV4	✓	✓	SC-TSV4	TSVC	v
	TSV 5.7	TSV4	✓	✓	SC-TSV5	TSVC	v
Astra	Astra EV 3.0	AEV30	✓		SC-AEV30	AEV30C	v
	Astra EV 3.6	AEV36	✓		SC-AEV36	AEV36C	v
	Astra EV 4.2	AEV42	✓		SC-AEV42	AEV42C	v
	Astra EV 4.8	AEV48	✓		SC-AEV48	AEV4854C	v
	Astra EV 5.4	AEV54	✓		SC-AEV54	AEV4854C	v
	Astra 3.5	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.0	ASR	✓	✓	SC-ASR	ASRC	v
	Astra 4.5	ASW	✓	✓	SC-ASW	ASWC	v
	Astra 5.0	ASW	✓	✓	SC-ASW	ASWC	v
Camlog	Camlog 3.3	CAM33	✓		SC-CAM33	CAMC	v
	Camlog 3.8	CAM38	✓		SC-CAM38	CAMC	v
	Camlog 4.3	CAM43	✓		SC-CAM43	CAMC	v
	Camlog 5.0	CAM50	✓		SC-CAM5060	CAM5060C	v
	Camlog 6.0	CAM60	✓		SC-CAM5060	CAM5060C	v
Magazen	AnyOne	MAO	✓	✓	SC-DER	MAOC	v
	Anyridge	MARR	✓	✓	SC-MARR	MARRC	v
Dentis	S-clean	DET	✓	✓	SC-DER	DERC	v
Dio	UF	SURD	✓	✓	SC-SURO	SURDC	v
Neo	IS	DERN	✓	✓	SC-SURO	DERNC	v
MIS	C1 3.3	C1M	✓		SC-C1M	C1MC	v
	C1 3.75 / 4.2	C1R	✓		SC-C1R	C1C	v
	C1 5.0	C1W	✓		SC-C1W	C1C	v

Geo Product List (Scanbody / Multibase / CMFit / 3D Analog / CCM Abutment / Sirona Ti-base)

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (Sea Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base							
101	Osstem	TS	3.0/3.5	0101-SUM	Single	SC-SUM	GMS-SUMS	LL2-SUM13-H2	LS2-SUM-H2	CMFIT10-SUM-H2	CMFIT14-SUM-H2	3D-SUM	AC-SUM5	SUM-H2	SRN2-SUM-H2							
					Bridge			LL2-SUM20-H2								SUM-N2	-					
							GMS-SUML	LL2-SUM13-N2	LS2-SUM-N2	CMFIT10-SUM-N2												
102			4.0/4.5/5.0- 0102-SURO	Single	SC-SURO	GMS-SUROS	LL2-SURO13-H2	LS2-SURO-H2	CMFIT10-SURO-H2	CMFIT14-SURO-H2	3D-SURO	AC-SURO5	SURO-H2	SRN2-SURO-H2								
				Bridge			LL2-SURO20-H2								SURO22-H2							
							LL2-SURO30-H2								SURO30-H2	SRN2-SURO26-H2						
103		US	3.5	0103-EXM	Single	SC-EXM	GMS-EXMS	LL2-EXM-H2	LS2-EXM-H2	CMFIT10-EXM-H2	CMFIT14-EXM-H2	3D-EXM	AC-EXM5	EXM-H2	-							
					Bridge			LL2-EXM-N2							EXM-N2	-						
									GMS-EXSL	LL2-EXR-H2				LS2-EXR-H2	CMFIT10-EXR-H2	EXR-H2	SRN2-EXR-H2					
104				4.1	0104-EXR	Single	SC-EXR	GMS-EXRL	LL2-EXR-N2		CMFIT10-EXR-N2	CMFIT14-EXR-H2	3D-EXR	AC-EXR5	EXR-N2							
						Bridge			GMS-EXWS	LL2-EXW-H2	LS2-EXW-H2				CMFIT10-EXW-H2	EXW-H2						
									GMS-EXWL	LL2-EXW-N2					CMFIT10-EXW-H2	EXW-N2	-					
105	4.8		0105-EXW	Single	SC-EXW	GMS-ISRS	LL2-ISR-O2	-	CMFIT10-ISR-O2	CMFIT14-ISR-O2	3D-ISR	AC-EXR5	ISR-O2	-								
				Bridge			GMS-ISRL	LL2-ISR-N2					ISR-N2	-								
							GMS-ISWS	LL2-ISW-O2					CMFIT10-ISW-O2	ISW-O2	-							
106			6.0	0106-ISR	Single	SC-ISRW	GMS-ISWL	LL2-ISW-N2		CMFIT10-ISW-O2	CMFIT14-ISW-O2	3D-ISW	AC-EXR5	ISW-N2	-							
					Bridge																	
107	3.7/4.0/4.5 (Universal)	0201-DER		Single	SC-DER	GMS-DERS	LL2-DER13-H2	LS2-DER-H2	CMFIT10-DER-H2	CMFIT14-DER-H2	3D-DER	AC-DER5	DER-H2	SRN2-DER-H2								
				Bridge			LL2-DER20-H2						DER20-H2	SRN2-DER26-H2								
							LL2-DER30-H2						DER30-H2	-								
201		Super-Line	3.7/4.0/4.5 (Universal)	0201-DER	Single	SC-DER	GMS-DERL	LL2-DER13-N2	LS2-DER-N2	CMFIT10-DER-N2	CMFIT14-DER-N2	3D-DER	AC-DER5	DER-N2	-							
					Bridge			LL2-DER20-N2						DER20-N2	-							
								LL2-DER30-N2						DER30-N2	-							
202	NR-Line		3.1/3.6/4.3 (Standard)	0202-DNR	Single	OS-DNR	-	-	CMFIT10-DNR-R	-	-	AC-DNR5	-	-								
301		AnyRidge	3.5/4.0/4.5 (Universal)	0301-MARR	Single	SC-MARR	GMS-MARRS	LL2-MARR-H2	LS2-MARR-H2	CMFIT10-MARR-H2	CMFIT14-MARR-H2	3D-MARR	AC-MARR5	MARR-H2	SRN2-MARR-H2							
					Bridge			LL2-MARR30-H2						MARR-N2	-							
								LL2-MARR20-N2														
302	AnyOne		3.5/4.0/4.5/5.0 (Universal)	0302-DERA	Single	SC-DER	GMS-DERS	LL2-DERA13-H2	LS2-DERA-H2	CMFIT10-DERA-H2	CMFIT14-DERA-H2	3D-DER	AC-MAO5	-	SRN2-DERA-H2							
					Bridge			LL2-DERA20-H2						CMFIT10-DERA-N2								
								LL2-DERA30-H2						CMFIT10-DERA-N2								
303		EZ-Plus	3.3	0303-MEPM	Single	SC-MEPM	-	-	CMFIT10-MEPM-E	-	-	AC-MEPM5	-	-								
					Bridge		LL2-MEPR-H2			CMFIT10-MEPR-E	-	AC-MEPR5	-	-								
							LL2-MEPR-N2			CMFIT10-MEPR-E	-	AC-MEPR5	-	-								
304	EX-feel		4.8 (Regular)	0305-ITRK	Single	SC-ITRK	-	-	CMFIT10-ITRK-O2	CMFIT14-ITRK-O2	3D-ITR	AC-EXR5	ITRK-O2	-								
					Bridge		LL2-ITRK-N2						CMFIT10-ITRK-N2		ITRK-N2	-						
305		Mini	3.0/3.25	0306-MMN	Single	SC-MMN	-	-	CMFIT10-MMN-H	-	-	AC-MMN5	-	-								
401	SM		3.8/4.1	0401-DSMM	Single	SC-DSMM	-	-	CMFIT10-DSMM-T	-	-	AC-DSMM5	-	-								
402		Dio	4.5/5.0	0402-DSMR	Single	SC-DSMR	-	-	CMFIT10-DSMR-T	-	-	AC-DSMR5	-	-								

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (See Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base
403	Dio	UF / UFII	3.0/3.3 (Narrow)	0403-DN30	Single	SC-MMN	-	-	-	CMFIT10-DN30-H	-	-	AC-DN305	-	-
404			3.8/4.0/4.5/5.0 (Standard)	0404-SURD	Single	SC-SURO	GMS-SUROS	LL2-SURD-H2 LL2-SURD20-H2 LL2-SURD30-H2 LL2-SURD-N2 LL2-SURD20-N2 LL2-SURD30-N2	-	CMFIT10-SURD-H2	CMFIT14-SURD-H2	3D-SURO	AC-SURD5	SURD-H2	-
					Bridge		GMS-SUROL	LL2-DERN-H2 LL2-DERN20-H2 LL2-DERN-N2 LL2-DERN20-N2		CMFIT10-SURD-N2	CMFIT14-SURD-N2			SURD-N2	-
501	NeoBiotech	IS	3.5/4.0/4.5/5.0-	0501-DERN	Single	SC-SURO	GMS-DERNS	LL2-DERN-H2 LL2-DERN20-H2 LL2-DERN-N2	-	CMFIT10-DERN-H2	CMFIT14-DERN-H2	3D-DER	AC-DERN5	DERN-H2	SRN2-DERN-H2
601	Dentis	S-Clean/OneQ	3.9/4.2/4.7/5.2		Bridge		GMS-DERNL	LL2-DERN20-N2		CMFIT10-DERN-N2	CMFIT14-DERN-N2			DERN-N2	SRN2-DERN20-H2
				Single	GMS-DERS	LL2-DERT13-H2 LL2-DERT13-N2	CMFIT10-DERT-H2	CMFIT14-DERT-H2	-	-					
700	Straumann	BoneLevel	2.9 (SC)	0700-BLS	Single	-	GMS-BLSL	LL2-BLS-R2	-	CMFIT10-BLS-R2	-	3D-BLS	AC-BLS5	-	-
701			3.3 (NC)	0701-BLM	Single	SC-BLM	GMS-BLMS	LL2-BLM-R2 LL2-BLM20-R2 LL2-BLM-H2	-	CMFIT10-BLM-R2	CMFIT14-BLM-R2	3D-BLM	AC-BLM5	BLM-R2	SRN2-BLM-R2
					Bridge		GMS-BLML	LL2-BLM20-N2		CMFIT10-BLM-R2	CMFIT14-BLM-R2			BLM-N2	SRN2-BLM-R2
702	Straumann	Tissue Level	4.1/4.8 (RC)	0702-BLR	Single	SC-BLR	GMS-BLRS	LL2-BLR-R2 LL2-BLR20-R2 LL2-BLR30-R2	-	CMFIT10-BLR-R2	CMFIT14-BLR-R2	3D-BLR	AC-BLR5	BLR-R2	SRN2-BLR-R2
703			3.5 (NN)	0703-IT35	Bridge		GMS-BLRL	LL2-BLR-N2 LL2-BLR20-N2		CMFIT10-BLR-R2	CMFIT14-BLR-R2			BLR-N2	SRN2-BLR-R2
					Single	SC-IT35	-	-	CMFIT10-IT35-O	-	-	AC-IT355	-	-	
704	Straumann	Tissue Level	4.8 (RN)	0704-ITR	Single	SC-ITRW	GMS-ITRS	LL2-ITR-O2 LL2-ITR20-O2	-	CMFIT10-ITR-O2	CMFIT14-ITR-O2	3D-ITR	AC-IT5	ITR-O2	SRN2-ITR-O2
705			6.5 (WN)	0705-ITW	Bridge		GMS-ITRL	LL2-ITR20-N2		CMFIT10-ITR-O2	CMFIT14-ITR-O2			ITR-N2	SRN2-ITR-O2
					Single	SC-ITRW	GMS-ITWS	LL2-ITW-O2 LL2-ITW20-O2	CMFIT10-ITW-O2	CMFIT14-ITW-O2	3D-ITW	SRN2-ITW-O2			
710	Nobel Biocare	BLX	3.75/4.5 (RB)	0710-BLX	Single	-	GMS-BLXL	LL2-BLX-T2	-	CMFIT10-BLX-T2	CMFIT14-BLX-T2	3D-BLX	AC-BLX5	-	-
711			5.5/6.5 (WB)	0711-BLXW	Single	-		CMFIT10-BLXW-T2	CMFIT14-BLXW-T2	-	-				
801	Nobel Biocare	Active	3.0	0801-AC30	Single	OS-AC30	-	LL2-AC30-T2 LL2-AC30-N2	LS2-AC30-T2	CMFIT10-AC30-T2	-	3D-AC30	AC-AC305	-	-
802			3.5 (NP)	0802-ACN	Single	SC-ACN	GMS-ACNS	LL2-ACN-H2 LL2-ACN20-H2 LL2-ACN-N2	-	CMFIT10-ACN-H2	CMFIT14-ACN-H2	3D-ACN	AC-ACN5	ACN-H2	SRN2-ACN-H2
					Bridge		GMS-ACNL	LL2-ACN20-N2		CMFIT10-ACN-H2	CMFIT14-ACN-H2			ACN-N2	SRN2-ACN-H2
803	Nobel Biocare	Active	4.3/5.0 (RP)	0803-ACR	Single	SC-ACR	GMS-ACRS	LL2-ACR-H2 LL2-ACR20-H2 LL2-ACR-N2	-	CMFIT10-ACR-H2	CMFIT14-ACR-H2	3D-ACR	AC-ACR5	ACR-H2	SRN2-ACR-H2
804			5.5 (WP)	0804-ACW	Bridge		GMS-ACRL	LL2-ACR20-N2		CMFIT10-ACR-H2	CMFIT14-ACR-H2			ACR-N2	SRN2-ACR-H2
					Single	SC-ACW	GMS-ACWS GMS-ACWL	LL2-ACW-H2 LL2-ACW-N2	CMFIT10-ACW-H2	CMFIT14-ACW-H2	3D-ACW	-	-		
805	Nobel Biocare	Branemark	3.5 (NP)	0805-EXMB	Single	SC-EXMB	-	LL2-EXMB-H2 LL2-EXMB-N2	LS2-EXMB-H2	CMFIT10-EXMB-H2	CMFIT14-EXMB-H2	3D-EXM	AC-EXMB5	EXMB-H2	-
806			4.1 (RP)	0806-EXRB	Bridge	SC-EXRB	-	LL2-EXRB-H2 LL2-EXRB-N2	-	CMFIT10-EXRB-H2	CMFIT14-EXRB-H2	3D-EXR	AC-EXRB5	EXRB-H2	SRN2-EXRB-H2
					Single			LL2-EXRB-N2		CMFIT10-EXRB-H2	CMFIT14-EXRB-H2			EXRB-N2	SRN2-EXRB-H2
807			5.1 (WP)	0807-EXWB	Single	SC-EXWB	-	LL2-EXWB-H2 LL2-EXWB-N2	LS2-EXWB-H2	CMFIT10-EXWB-H2	CMFIT14-EXWB-H2	3D-EXW	AC-EXWB5	EXWB-H2	-

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMF fit 10mm	CMFIT 14mm	3D Analog	Screw (See Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base		
808	Nobel Biocare	Replace Select	3.5 (NP)	0808-RPM	Single	SC-RPM	GMS-RPMS	LL2-RPM-E2	LS2-RPM-E2	CMFIT10-RPM-E2	CMFIT14-RPM-E2	3D-RPM	AC-RPM5	RPM-E2	SRN2-RPM-E2		
					Bridge		GMS-RPML	LL2-RPM-N2				PRM-N2					
809			4.3 (RP)	0809-RPR	Single	SC-RPR	GMS-RPRS	LL2-RPR-E2	LS2-RPR-E2	CMFIT10-RPR-E2	CMFIT14-RPR-E2	3D-RPR	AC-RPR5	RPR-E2	SRN2-RPR-E2		
					Bridge		GMS-RPRL	LL2-RPR-N2				RPR-N2					
810			5 (WP)	0810-RPW	Single	SC-RPW	GMS-RPWS	LL2-RPW-E2	LS2-RPW-E2	CMFIT10-RPW-E2	CMFIT14-RPW-E2	3D-RPW		RPW-E2	SRN2-RPW-E2		
					Bridge		GMS-RPWL	LL2-RPW-N2				RPW-N2					
811	Dentsply	Xive	6.0	0811-RP60	Single	SC-RP60	GMS-RP60S	LL2-RP60-E2	LS2-RP60-E2	CMFIT10-RP60-E2	CMFIT14-RP60-E2	3D-RP60	AC-XV35	-	-		
					Bridge		GMS-RP60L	LL2-RP60-N2				-					
901			3.0	0901-XV3	Single	SC-XV3	GMS-XV3L	LL2-XV3-H2	LS2-XV3-H2	CMFIT10-XV3-H2	-	3D-XV3	-	-			
					Bridge		GMS-XV3S	LL2-XV3-N2				-					
902			3.4	0902-XVM	Single	SC-XVM	GMS-XVMS	LL2-XVM-H2	LS2-XVM-H2	CMFIT10-XVM-H2	CMFIT14-XVM-H2	3D-XVM	AC-XV5	XVM-H2	SRN2-XVM-H2		
					Bridge		GMS-XVML	LL2-XVM-N2				XVM-N2					
903			3.8	0903-XVR	Single	SC-XVR	GMS-XVRS	LL2-XVR-H2	LS2-XVR-H2	CMFIT10-XVR-H2	CMFIT14-XVR-H2	3D-XVR		XVR-H2	SRN2-XVR-H2		
					Bridge		GMS-XVRL	LL2-XVR-N2				XVR-N2					
904			4.5	0904-XVW	Single	SC-XVW	GMS-XVWS	LL2-XVW-H2	LS2-XVW-H2	CMFIT10-XVW-H2	CMFIT14-XVW-H2	3D-XVW	AC-XV5	XVW-H2	SRN2-XVW-H2		
					Bridge		GMS-XVWL	LL2-XVW-N2				XVW-N2					
905	5.5	0905-XV55	Single	SC-XVW	GMS-XVWS	LL2-XV55-H2	LS2-XV55-H2	CMFIT10-XV55-H2	CMFIT14-XV55-H2	3D-XV5	-	-					
			Bridge		GMS-XVWL	LL2-XV55-N2				-							
906	Ankylos	3.4/4.5/5.5 (Universal)	0906-AKLL	Single	OS-AKLL	-	LL2-AKLL-T		CMFIT10-AKLL-T	CMFIT14-AKLL-T	3D-AKLL	-	-	SRN2-AKLL 1-T			
				Bridge			LL2-AKLL20-T				-		SRN2-AKLL 2-T				
			DNF3-N	LL2-AKLL-N	CMFIT10-AKLL-N		-		SRN2-AKLL 1-N								
				LL2-AKLL20-N					SRN2-AKLL 2-N								
1001	Biomet 3i	Certain	3.4	1001-3ICM	Single	SC-3ICM	GMS-3ICMS	LL2-3ICM-H2	LS2-3ICM-H2	CMFIT10-3ICM-H2	-	3D-3ICM	AC-3IC5	3ICM-H2	-		
					Bridge		GMS-3ICML	LL2-3ICM-N2				3ICM-N2					
1002			4.1	1002-3ICR	Single		GMS-3ICRS	LL2-3ICR-H2	LS2-3ICR-H2	CMFIT10-3ICR-H2	CMFIT14-3ICR-H2	3D-3ICR		3ICR-H2	-		
					Bridge		GMS-3ICRL	LL2-3ICR-N2				3ICR-N2					
1003			5.0	1003-3ICW	Single	SC-3ICR	GMS-3ICRS	LL2-3ICW-H2	LS2-3ICW-H2	CMFIT10-3ICW-H2	CMFIT14-3ICW-H2	3D-3ICW	-	-			
					Bridge		GMS-3ICRL	LL2-3ICW-N2				-					
1004			6.0	1004-3IC60	Single		GMS-3ICRS	LL2-3IC60-H2	LS2-3IC60-H2	CMFIT10-3IC60-H2	CMFIT14-3IC60-H2	3D-3IC60	-	-			
					Bridge		GMS-3ICRL	LL2-3IC60-N2				-					
1005	External / Osseotite		3.4	1005-EX3M	Single	SC-EX3M	-	LL2-EX3M-H2	-	CMFIT10-EX3M-H2	-	-	AC-EXR5	-	-		
					Bridge			LL2-EX3M-N2				-		-			
1006			4.1	1006-EXR	Single		GMS-EXRS	LL2-EXR-H2	LS2-EXR-H2	CMFIT10-EXR-H2	CMFIT14-EXR-H2	3D-EXR		EXR-H2	SRN2-EXR-H2		
					Bridge		GMS-EXRL	LL2-EXR-N2		CMFIT10-EXR-N2		EXR-N2					
1007			5.0	1007-EXT	Single	SC-EXR	GMS-EXRS	LL2-EXT-H2	LS2-EXT-H2	CMFIT10-EXT-H2	CMFIT14-EXT-H2	-	EXT-H2	-			
					Bridge		GMS-EXRL	LL2-EXT-N2				EXT-N2					
1008			6.0	1008-EXTW	Single		GMS-EXRS	LL2-EXTW-H2	-	CMFIT10-EXTW-H2	CMFIT14-EXTW-H2	-	-	-	-		
					Bridge		GMS-EXRL	LL2-EXTW-N2				-					
1101	Implant Direct	Legacy	3.0	1101-TSMV	Single	OS-TSMV	-	LL2-TSMV-H2	-	CMFIT10-TSMV-H2	-	3D-TSMV	AC-TSMV5	-	-		
1201	Zimmer	Tapered Screw-Vent	3.5	1201-TSV3	Single	SC-TSV3	GMS-TSV3S	LL2-TSV3-H2		CMFIT10-TSV3-H2	CMFIT14-TSV3-H2	3D-TSV3	AC-TSV5	TSV3-H2	SRN2-TSV3-H2		
					Bridge			LL2-TSV320-H2				TSV3-N2				TSV3-N2	
							GMS-TSV3L	LL2-TSV320-N2				TSV3-H2					
1202			4.5	1202-TSV4	Single	SC-TSV4	GMS-TSV4S	LL2-TSV4-H2	LS2-TSV4-H2	CMFIT10-TSV4-H2	CMFIT14-TSV4-H2	3D-TSV4	TSV4-H2	SRN2-TSV4-H2			
					Bridge		GMS-TSV4L	LL2-TSV420-H2				TSV4-N2			TSV4-N2		
							GMS-TSV4L	LL2-TSV420-N2				-					
1203			5.7	1203-TSV5	Single	SC-TSV5	GMS-TSV5S	LL2-TSV5-H2	LS2-TSV5-H2	CMFIT10-TSV5-H2	CMFIT14-TSV5-H2	3D-TSV5	TSV5-H2	SRN2-TSV5-H2			
	Bridge	GMS-TSV5L			LL2-TSV5-N2					TSV5-N2			TSV5-N2				
1204	Spline	3.7/4.0			1204-SP40		Single	SC-SP40	-	LL2-SP40-H	-	-	CMFIT14-SP40-H	-	SP40-H	-	
1205			5.0	1205-SP50		Single	-		LL2-SP50-H	-	-	CMFIT14-SP50-H	-	-	-		

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFit 10mm	CMFit 14mm	3D Analog	Screw (Sea Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base		
1301	Osseospeed		3.0 (Yellow)	1301-AS3	Single	SC-AS3	- GMS-AS3L	LL2-AS3-T2	-	CMFIT10-AS3-T2	-	3D-AS3	AC-AS35	AS3-T2	-		
				1302-ASR	Single	SC-ASR	GMS-ASRS	LL2-ASR-H2	LS2-ASR-H2	CMFIT10-ASR-H2	CMFIT14-ASR-H2	3D-ASR	AC-ASR5	ASR-H2	SRN2-ASR-H2		
1302			Bridge		LL2-ASR20-H2			LL2-ASR20-N2									
			Single		LL2-ASW-H2									LS2-ASW-H2		CMFIT10-ASW-H2	CMFIT14-ASW-H2
1303			4.5/5.0 (Lilac)	1303-ASW	Bridge	LL2-ASW20-H2	LL2-ASW-N2	ASW-N2									
1304					3 (Green)	1304-AEV30	Single Bridge	GMS-AEV30S GMS-AEV30L	LL2-AEV30-T2 LL2-AEV30-N2	-	CMFIT10-AEV30-T2	-	3D-AEV30		AC-AEV305		
1305	Astra Tech	EV	3.6 (Purple)	1305-AEV36	Single	SC-AEV36	GMS-AEV36S	LL2-AEV36-T2	-	CMFIT10-AEV36-T2	CMFIT14-AEV36-T2	3D-AEV36	AC-AEV365	-	-		
					Bridge			LL2-AEV36-N2						LL2-AEV3620-N2			
1306			4.2 (Yellow)	1306-AEV42	Single	SC-AEV42	GMS-AEV42S	LL2-AEV42-T2	-	CMFIT10-AEV42-T2	CMFIT14-AEV42-T2	3D-AEV42	AC-AEV425	-	-		
					Bridge			LL2-AEV42-N2						LL2-AEV4220-N2			
1307			4.8 (Blue)	1307-AEV48	Single	SC-AEV48	GMS-AEV48S	LL2-AEV48-T2	-	CMFIT10-AEV48-T2	CMFIT14-AEV48-T2	3D-AEV48	AC-AEV48545	-	-		
					Bridge			LL2-AEV48-N2						LL2-AEV4820-N2			
1308	Camlog		5.4 (Brown)	1308-AEV54	Single	SC-AEV54	GMS-AEV54S	LL2-AEV54-T2	-	CMFIT10-AEV54-T2	CMFIT14-AEV54-T2	3D-AEV54	AC-AEV545	-	-		
					Bridge			LL2-AEV54-N2						LL2-AEV5420-N2			
1401			Camlog		3.3	1401-CAM33	Single Bridge	SC-CAM33	-	LL2-CAM33-T2	-	CMFIT10-CAM33-T2	-	3D-CAM33	AC-CAM5	-	-
1402							Single Bridge			LL2-CAM33-N2						LL2-CAM38-T2	CMFIT10-CAM38-T2
1403					4.3	1403-CAM43	Single Bridge	SC-CAM43	-	LL2-CAM43-T2	-	CMFIT10-CAM43-T2	CMFIT14-CAM43-T2	3D-CAM43	AC-CAM50605	-	-
1404							Single			LL2-CAM43-N2						LL2-CAM50-T2	CMFIT10-CAM50-T2
1405	6.0	1405-CAM60			Single Bridge	SC-CAM5060	-	LL2-CAM60-T2	-	CMFIT10-CAM60-T2	CMFIT14-CAM60-T2	3D-CAM60	AC-CAM50605	-	-		
1406					3.3			1406-CN33						Single	-	GMS-CN33L	-
1407	Conelog		3.8/4.3	1407-CN3843	Single	-	GMS-CN3843L	-	-	CMFH14-CN3843-E	3D-CN3843	-	AC-CN5	-	-		
1408					5.0	1408-CN50	Single	-	GMS-CN50L	-	-	CMFH14-CN50-E	3D-CN50	AC-CN505	-	-	
1501	MIS	Seven	3.3 (Narrow)	1501-MS33	Single Bridge	SC-MS33	-	LL2-MS33-H	-	CMFIT10-MS33-H	-	3D-MS33	AC-MS335	-	-		
1502					Single Bridge			LL2-MS33-N						LL2-MSR-H	CMFIT10-MSR-H	CMFIT14-MSR-H	3D-TSV3
1503			5.0/6.0 (Wide)	1503-MSW	Single Bridge	SC-TSV4	GMS-TSV3L GMS-TSV4S	LL2-MSR-N	-	CMFIT10-MSW-H	CMFIT14-MSW-H	3D-TSV4	AC-MS5	-	-		
1504					Single Bridge			LL2-MSW-N						LL2-MSW-H	CMFIT10-MSW-H	CMFIT14-MSW-H	3D-TSV4
1505			C1		3.3	1504-C1M	Single Bridge	SC-C1M	GMS-C1MS	LL2-C1M-T	-	CMFIT10-C1M-T	CMFIT14-C1M-T	3D-C1M	AC-C1M5	C1M-T	-
1506							Single Bridge			LL2-C1M-N						LL2-C1R-T	CMFIT10-C1R-T

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMFIt 10mm	CMFIt 14mm	3D Analog	Screw (See Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base
1601	SIC	ace / SIC max	3.3	1601-SICM	Single	SC-SIC	-	LL2-SICM-H	-	CMFIT10-SICM-H	CMFIT14-SICM-H	3D-SICM	AC-SIC5	SICM-H	-
1602				1602-SICR	Bridge			LL2-SICM-N		CMFIT10-SICR-H	CMFIT14-SICR-H	3D-SICR		SICR-H	-
1701	Thommen	SPI	3.5	1701-SPI35	Single	SC-SPI35	-	LL2-SPI35-H	-	CMFIT10-SPI35-H	-	3D-SPI35	AC-SPI355	-	-
1702				1702-SPI40	Bridge			LL2-SPI35-N		CMFIT10-SPI40-H	CMFIT14-SPI40-H	3D-SPI40		-	-
1703			4.5	1703-SPI45	Single	SC-SPI45	-	LL2-SPI45-H	-	CMFIT10-SPI45-H	CMFIT14-SPI45-H	3D-SPI45	AC-SPI5	-	-
1704				1704-SPI50	Bridge			LL2-SPI45-N		CMFIT10-SPI50-H	CMFIT14-SPI50-H	3D-SPI50		-	-
1705			6.0	1705-SPI60	Single			LL2-SPI60-H	-	CMFIT10-SPI60-H	CMFIT14-SPI60-H	-		-	-
1801			3.7	1801-BEGO37	Single	SC-BEGO	-	LL2-BEGO37-H	-	CMFIT10-BEGO37-H	-	3D-BEGO37	AC-BEGO5	-	-
1802	Bego	Bego	4.1	1802-BEGO41	Bridge			LL2-BEGO37-N		CMFIT10-BEGO37-N	-	3D-BEGO41		-	-
1803				1803-BEGO45	Single			LL2-BEGO41-H	-	CMFIT10-BEGO41-H	-	3D-BEGO45		-	-
1804			5.5	1804-BEGO55	Single			LL2-BEGO41-N		CMFIT10-BEGO45-H	CMFIT14-BEGO45-H	3D-BEGO45		-	-
1901	Mytis	Arrow	3.7/4.0	1901-AR3740	Single	-	GMS-AR3740L	-	-	CMFIT10-BEGO55-H	CMFIT14-BEGO55-H	3D-BEGO55	AC-AR5	-	-
1902			4.7/6.0	1902-AR4760	Single	OS-AR4760				CMFIT10-AR3740-H	-	-		-	-
2001	KYOCERA	POI	3.7	2001-POI37	Single	SC-POI37	-	LL2-POI37-H	-	CMFIT10-POI37-H	-	-	AC-POI375	-	-
2002			4.2	2002-POI42	Single	SC-POI42	-	LL2-POI42-H	-	CMFIT10-POI42-H	-	-		-	-
2101	Biohorizon	Laser Lok	3.0	2101-BH30	Single	OS-BH30	-	LL2-BH30-H2	-	CMFIT10-BH30-H2	-	3D-BH30	AC-BH305	-	-
2201	Anthogyr	AXIOM Bone Level	3.4/4.0/4.6/5.2 (Universal)	2201-ANG	Single	SC-ANG	GMS-ANGS	LL2-ANG-E	-	CMFIT10-ANG-E	-	3D-ANG	AC-ANG5	ANG-E	-
2301	WARANTEC	Oneplant	4.3/4.8/5.3 (Standard)	2301-WOP	Single	SC-WOP	-	LL2-ANG-N	-	CMFIT10-WOP-O	-	-	AC-WOP5	-	-
2302		IU	3.3/3.6 (Mini)	2302-SUMW	Single	-	GMS-SUMWL	LL2-SUMW-H2	-	CMFIT10-SUMW-H2	-	3D-SUMW	AC-SUMW5	-	-
2401	Tekta	In-Kone	3.5/4.0/4.5/5.5 (Standard)	2401-TEK	Single	SC-TEK	-	-	-	CMFIT10-TEK-H	-	3D-TEK	AC-TEK5	TEK-H	-
2501	Keystone	PrimaComex / Max-Trilobe	3.5	2501-KEP35	Single	SC-KEP35	GMS-KEP35S	LL2-KEP35-T2	LS2-KEP35-T2	CMFIT10-KEP35-T2	CMFIT14-KEP35-T2	3D-KEP35	AC-KEP5	-	-
2502				2502-KEP41	Bridge	SC-KEP41	GMS-KEP35L	LL2-KEP35-N2		CMFIT10-KEP41-T2	CMFIT14-KEP41-T2	3D-KEP41		-	-
2503			5.0/6.0	2503-KEP50	Single	SC-KEP50	GMS-KEP41S	LL2-KEP41-T2	LS2-KEP41-T2	CMFIT10-KEP41-T2	CMFIT14-KEP41-T2	3D-KEP50	AC-KEP5	-	-
2504				2504-KEX50	Bridge		GMS-KEP41L	LL2-KEP41-N2	LS2-KEP50-T2	CMFIT10-KEP50-T2	CMFIT14-KEP50-T2	3D-KEX50		-	-
2505			8.0	2505-KEX65	Single	SC-KEP50	GMS-KEP50S	LL2-KEP50-T2	LS2-KEX50-T2	-	CMFIT14-KEX50-T2	3D-KEX50	-	-	-
2601	Neodent	GM	3.5/4.3/5.0 (Standard)	2601-NCM	Single	SC-NCM	-	LL2-KEX50-N2	LS2-KEX65-T2	CMFIT14-KEX65-T2	-	3D-NCM	-	-	-
2602				2602-NCG	Single	SC-NCG	GMS-NCGL	LL2-NCG-H2	-	CMFIT10-NCG-H2	CMFIT14-NCG-H2	3D-NCG	AC-NCG5	-	-

NO	Company	Implant system	Platform	Library Code	Type	SC(OS) Scanbody	GMS Scanbody	LL-Multibase	LS-Multibase	CMF fit 10mm	CMF fit 14mm	3D Analog	Screw (sea Pack)	CCM Abutment (Co-Cr base)	Sirona Ti-base
2701	BICON	BICON Implant	2.0	2701-BC20	Bridge	SC-BC20	-	LL2-BC20-N	-	CMFIT10-BC20-N	-	-	-	-	-
2702			2.5	2702-BC25	Bridge	SC-BC25	-	LL2-BC25-N	-	CMFIT10-BC25-N	-	-	-	-	-
2703			3.0	2703-BC30	Bridge	SC-BC30	-	LL2-BC30-N	-	CMFIT10-BC30-N	-	-	-	-	-
2801	Medentika	IPS	3.5/4.0/4.5/5.0	2801-MTR	Single	SC-MTR	-	-	-	CMFIT10-MTR-R	-	-	AC-MTR5	-	-
2901	GC	Aadva Implant	4.0/5.0 (Regular / Wide)	2901-GAD	Single	SC-GAD	-	-	-	CMFIT10-GAD-H	-	-	AC-GAD5	-	-
3001	BIOTECH	KONTACT	3.0	3001-BK30	Single Bridge	SC-BK30	-	LL2-BK30-T LL2-BK30-N	-	CMFIT10-BK30-T	-	3D-BK30	AC-BK305	-	-
3002			3.6/4.2/4.8/5.4	3002-BK40	Single Bridge	SC-BK40	-	LL2-BK40-T LL2-BK40-N	-	CMFIT10-BK40-T	-	3D-BK40	AC-BK405	-	-
3101	Medentis medical	ICX	3.75/4.1/4.8	3101-ICX	Single Bridge	SC-ICX	-	LL2-ICX-H LL2-ICX-N	-	CMFIT10-ICX-H	CMFIT14-ICX-H	3D-ICX	AC-ICX5	-	-
3201	B&B Dental	Dura-Vit	4.0/3.5/5.0	3201-BNB	Single	-	GMS-BNBL	LL2-BNB-H	-	CMFIT10-BNB-H	-	3D-BNB	AC-BNB5	-	-
3301	Bredent Medical	Sky	3.5 / 4.0	3301-SKR	Single	-	GMS-SKRL	LL2-SKR-T	-	CMFIT10-SKR-T	CMFIT14-SKR-T	3D-SKR	AC-SKR5	-	-
3401	Necoss	ProActive	3.5/4.0/4.5/5.0/5.5	3401-NEP	Single	-	GMS-NEPL	LL2-NEP-T	-	-	CMFIT14-NEP-T	3D-NEP	AC-NEP5	-	-
3501	Champions	(R)evolution	3.5/4.0/4.5/5.5	3501-CHR	Single	-	GMS-CHRL	-	-	-	CMFIT14-CHR-H	3D-CHR	AC-CHR5	-	-
3601	KYOCERA	FINESIA Bone Level	3.2 / 3.4	3601-FBLN	Single	-	GMS-FBLNL	-	-	CMFIT10-FBLN-H	-	3D-FBLN	AC-FBLN5	-	-
3602			3.7 / 4.2	3602-FBLR	Single	-	GMS-FBLRL	LL2-FBLR-H	-	CMFIT10-FBLR-H	-	3D-FBLR	AC-FBLR5	-	-

Geo MultiUnit List (Scanbody / Cylinder / 3D Analog)

NO	Compatible to Implant Systems	Platform	Library Code	Type	SC(OS) Scanbody	LB2-short	LB2-long	Temporary Cylinder	3D Analog	Screw(Sea)	CCM Abutment (Co-Cr base)
203	Dentium Screw Abutment	4.5	0203-MUDE45	Bridge	SC-MUDE	LB2-MUDE45	-	-	-	AC-MUDE5	-
204		5.5	0204-MUDE55			LB2-MUDE55					
706	Straumann TissueLevel	4.8	0706-MUIT	Bridge	SC-MUIT	LB2-MUIT	-	-	-	AC-MUIT5	-
707	Straumann Bonelevel SRA	3.4	0707-MUSR35	Bridge	SC-MUSR35	LB2-MUSR35	-	-	3D-MUSR35 3D-MUSR46	AC-MUSR5	-
708		4.6	0708-MUSR46		SC-MUSR46	LB2-MUSR46					
812	Nobel Multi Unit	4.8	0812-MUTB	Bridge	SC-MUTB	LB2-MUTB2	LB2-MUTBL 2	TA-MUTB	3D-MUTB	AC-MUTB5	MUTB-N2
907	Dentsply Ankylos Multi-Unit	mini	0907-AKLM	Bridge	SCB-16	LB2-AKLM	-	-	-	AC-AKL5	-
908		regular	0908-AKLR	Bridge		LB2-AKLR					
1206	Zimmer Tapered Screw-Vent Tapered Coping	4.5	1206-MUTSV	Bridge	SC-MUTSV	LB2-MUTSV	-	-	3D-MUTSV	AC-MUTSV5	-
108	Ossitem Multi Abutment	4.8	0108-MUTO	Bridge	SC2-MUTO	LB2-MUTO	LB2-MUTOL	-	-	AC-MUTO5	-

Geo Transfer Scanbody List

NO	Compatible to Implant Systems	Platform	Type	Section				
101	Osstem Implant TS Mini	3.0/3.5	Single	Diameter GH / Height	4.5mm	5.0mm		
				2.0 / 5.0	GTS-SUM4020	GTS-SUM4520		
				4.0 / 7.0	GTS-SUM4040	GTS-SUM4540		
				6.0 / 9.0	GTS-SUM4060	GTS-SUM4560		
102	Osstem Implant TS Standard	4.0/4.5/5.0-	Single	Diameter GH / Height	4.5mm	5.0mm	6.0mm	7.0mm
				2.0 / 5.0	GTS-SURO4520	GTS-SURO5020	GTS-SURO6020	GTS-SURO7020
				4.0 / 7.0	GTS-SURO4540	GTS-SURO5040	GTS-SURO6040	GTS-SURO7040
				6.0 / 9.0	GTS-SURO4560	GTS-SURO5060	GTS-SURO6060	GTS-SURO7060
201	Dentium Super-Line Implantium	3.7/4.0/4.5-	Single	Diameter GH / Height	4.0mm	4.5mm	5.5mm	6.5mm
				2.0 / 5.0	GTS-DER4020	GTS-DER4520	GTS-DER5520	GTS-DER6520
				4.0 / 7.0	GTS-DER4040	GTS-DER4540	GTS-DER5540	GTS-DER6540
				6.0 / 9.0	GTS-DER4060	GTS-DER4560	GTS-DER5560	GTS-DER6560

Guide Tool	Compatible to Implant Systems	Diameter						
		4.0mm	4.5mm	5.0mm	5.5mm	6.0mm	6.5mm	7.0mm
	Osstem Implant TS Mini	GTST-S40-120H	GTST-S45-120H	-	-	-	-	-
	Osstem Implant TS Standard	-	GTST-S45-120H	GTST-S50-120H	-	GTST-S60-120H	-	GTST-S70-120H
	Dentium Super-Line Implantium	GTST-D40-127H	GTST-D45-127H	-	GTST-D55-127H	-	GTST-D65-127H	-

GeOMEDI



Digital dentistry has revolutionized the dental industry forever.
To provide a better CAD-CAM experience to our customers, we founded GeoMedi in 2009.

GeoMedi's core value is diversity, expertise, and customer satisfaction.

We have one of the most diverse compatible implant systems which meet most of your requirements.
We strive to develop more implant systems with superior quality.

(주) 지오메디

- ▣ 경기도 의왕시 이미로 40 (포일동, 인덕원IT밸리) D동 B112, 702~705, 709-1, 810-2, 811호
- ▣ Tel : 02-588-2871
- ▣ Fax : 02-6280-2849
- ▣ E-mail : sales@geomedi.com
- ▣ WebSite : www.geomedi.com

GeoMedi Co., Ltd.

- ▣ #B112, 702~705, 709-1, 810-2, 811, Unit D, Indeokwon IT Valley, 40, Imi-ro, Uiwang-si, Gyeonggi-do, Korea
- ▣ Tel : 82-2-588-2871
- ▣ Fax : 82-2-6280-2849
- ▣ E-mail : sales@geomedi.com
- ▣ WebSite : www.geomedi.com



GeoMedi

PREMIUM CAD/CAM PRODUCTS